



EXPERTS IN PHYSIOTHERAPY



2023

ALL THERAPY TABLES
CUSTOMIZED AND
MADE IN GERMANY

MADE
IN
GERMANY



WELCOME TO PINO

The leading supplier of professional physiotherapy products in Germany, Austria, Switzerland and the Netherlands.

PINO has been the reliable partner of physiotherapists and sports doctors in Germany, Austria, Switzerland and the Netherlands for almost 120 years. This enormous experience has resulted in products of exceptional quality that meet the highest demands of professional users. PINO always aims to drive new developments - developed by physiotherapists for physiotherapists! Our exclusive own brands offer you everything that the modern practice needs, from kinesiological tapes, massage products and training equipment to therapy beds.

Discover our range and let us convince you of the combination of tradition and modernity.

COMPANY HISTORY

Travel with us through our company chronology and discover the most important stages of our history.



Start of the production of massage oils.

1937



Foundation in Freudenstadt in the Black Forest.

1904



Relocation of the company headquarters to Hamburg.

1994



Start of own Production of therapy tables in Germany.

2010



Training of the 20,000 Therapist in the basic course PINOTAPE 1.

2014



Sales launch of PINOTAPE kinesiology tapes.

2006



Presentation of PINOFIT Physioboost massage gun.

2019

Check out all our products

THERAPY & TRAINING



GYMNASTIC BANDS IN PROFESSIONAL QUALITY

Various types & materials for every field of application.



	Gymnastic bands	Stretch bands
Material	Pure natural latex <ul style="list-style-type: none"> Made from pure natural latex (99% free from allergens) 	Skin-friendly textile material <ul style="list-style-type: none"> Washable at 30 °C
Resistances	Five different resistances * Black only available in 45 m length	Three different resistances
Lengths	 2 m 45 m 33 cm (loop)	 1 m 2.70 m 33 cm (loop)
Your benefits	<ul style="list-style-type: none"> 45 m bands: can be trimmed to size 2 m bands + loops: individually packed, suitable for resale 	<ul style="list-style-type: none"> Numbered loops: makes explaining & performing exercises easy Individually packed, suitable for resale

GYMNASTIC BAND MADE FROM NATURAL LATEX

Professional quality at an attractive price.



✓ MEDICAL PRODUCT CE

PINOFIT® gymnastic band Display

Includes 24 packs of PINOFIT® gymnastic band, 2 m each.

Classic

Item no. 45075

24 bands in classic colours, 8 red and 8 green bands, 4 yellow and 4 blue bands

Fresh

Item no. 44975

24 bands in fresh colours, 8 coral and 8 lime bands, 4 yellow and 4 azure bands



✓ MEDICAL PRODUCT CE

PINOFIT® Gymnastic band 2 m pack

- 2 m long, 14 cm wide
- Perfect for therapy and home use
- No trimming necessary
- Made from pure natural latex (99% free from allergenic latex proteins)
- Verified quality: more than 10,000 stretch repetitions without faults

Yellow

Item no. 45070 Resistance light, 2 m long

Red

Item no. 45071 Resistance medium, 2 m long

Lime

Item no. 45074 Resistance strong, 2 m long

Azure

Item no. 45078 Resistance extra strong, 2 m long

Coral

Item no. 44077 Resistance medium, 2 m long

Green


Item no. 45072 Resistance strong, 2 m long

Blue

Item no. 45073 Resistance extra strong, 2 m long

All colours freely combinable

✓ MEDICAL PRODUCT CE



PINOFIT® Gymnastic band 45 m roll

- Offers optimal stretching properties combined with linear resistance development
- All important muscle groups can be targeted and exercised effectively with the PINOFIT® gymnastic band
- For therapeutic training and rehabilitation of muscle and joint functions
- Made from pure natural latex (99% free from allergenic latex proteins)
- Without additives, therefore partly transparent
- Verified quality: more than 10,000 stretch repetitions without faults
- Professional quality
- Extra rip-proof and durable
- No adhesion despite the lack of powder

Yellow

Item no. 45001

Resistance light, 45 m long

Coral

Item no. 44007

Resistance medium, 45 m long

Red

Item no. 45002

Resistance medium, 45 m long

Lime

Item no. 45014

Resistance strong, 45 m long

Green

Item no. 45003

Resistance strong, 45 m long

Azure

Item no. 44008

Resistance extra strong, 45 m long

Blue

Item no. 45004

Resistance extra strong, 45 m long

Black

Item no. 45005

Resistance super strong, 45 m long



✓ MEDICAL PRODUCT CE

PINOFIT® Miniband

- Small gymnastics band for great effects
- The wider quality withstands highest pressure without any problems in comparison to traditional minibands.
- Small, looped gymnastic band
- Versatile for strength and stretch exercises
- Ideal for on the go

Yellow

Item no. 44650

Resistance light, approx. 33 x 6 cm

Lime

Item no. 44654

Resistance strong, approx. 33 x 6 cm

Blue

Item no. 44653

Resistance extra strong, approx. 33 x 6 cm

Red

Item no. 44651

Resistance medium, approx. 33 x 6 cm

Green

Item no. 44652



Resistance strong, approx. 33 x 6 cm

All colours freely combinable



STRETCH BAND MADE FROM TEXTILE MATERIAL

For static and dynamic stretching.



PINOFIT® Stretch Band

- For therapy, rehab & yoga
- Durable & long-lasting
- Washable at 30 °C
- 39% polyamide, 36% spandex and 25% polyester
- 8 handles with numbers for orientation
- Latex-free

Yellow


Item no. 44643 Resistance light, approx. 2.5 cm wide x 100 cm long

Coral

Item no. 44644 Resistance medium, approx. 4 cm wide x 100 cm long

Lime

Item no. 44645 Resistance strong, approx. 5 cm wide x 100 cm long



PINOFIT® Stretch Band XL

- Extra-long stretch band for even more possible uses
- 39% polyamide, 36% spandex and 25% polyester
- Washable at 30 °C
- 16 handles with numbers for orientation
- Latex-free

Yellow

Item no. 44656 Resistance light, approx. 2.5 cm wide x 270 cm long

Coral

Item no. 44648 Resistance medium, approx. 4 cm wide x 270 cm long

Lime

Item no. 44659 Resistance strong, approx. 5 cm wide x 270 cm long



Pack of 2

PINOFIT® Stretch Miniband

- Small, looped stretch band
- Ideal for exercises for activating and stabilising muscles
- 39% polyamide, 36% spandex and 25% polyester
- Washable at 30 °C
- Latex-free

Yellow

Item no. 44629 Resistance light, approx. 4 cm wide x 33 cm long, pack of 2

Coral

Item no. 44649 Resistance medium, approx. 4 cm wide x 33 cm long, pack of 2

Lime

Item no. 44669 Resistance strong, approx. 4 cm wide x 33 cm long, pack of 2



Our highlights for your training









- ⊕ Very comfortable thanks to skin-friendly textile material
- ⊕ Washable
- ⊕ Durable







OUR WAVE PRODUCTS

PINO developed the PINOFIT® Wave together with physiotherapists in 2013. The innovative wave profile of the fascia roller enables intensive effects during massage, supporting the release of fascial tension and adhesions. We now offer the perfect roller for every purpose and taste: choose your favourite(s) from a range of sizes, shapes, colours and hardness levels!



				
	Fascia Wave Roller	fascia roller Wave Pro	fascia roller Wave Stick	fascia roller Wave Mini
	Extra long & extra light	With ergonomic indentation	Handles allow individual application of pressure	Extra compact & handy
Hardness	●●○ Medium	●●○ Medium	●●● Hard	●●○ Medium
Available colours				
Size	45 cm long, 12 cm Ø	36.5 cm long, 13 cm Ø	Without handles approx. 28 cm long, ø approx. 8 cm	Approx. 16 cm long, ø approx. 6 cm
Application areas	Universal for any fascia training	Especially for fascia training on the back	Especially for fascia training while sitting or standing (e.g. in rehab exercise)	Especially for fascia training on soles of feet, calves & arms



		
	Fascia ball	Fascia duoball
	Multifunctional & compact	With ergonomic notching
Hardness	●●● Hard	●●● hard
Available colours		
Size	ø approx. 9.8 cm	Length approx. 15 cm
Application areas	<ul style="list-style-type: none"> • For targeted application on trigger points • Spot application 	<ul style="list-style-type: none"> • For targeted application in the neck area • Spot application

FASCIAL ROLLERS WITH WAVE PROFILE

With a wave for a more intensive training effect.



Our highlights for your application

- ⊕ Ergonomically shaped contour with perfectly shaped indentation to provide an ideal relief for the spine (and other parts of the body) during the workout
- ⊕ In tried-and-tested Wave quality
- ⊕ Perfect for back exercises
- ⊕ Fascia roller workout videos at: www.youtube.de/PINOFIT

TOP SELLER



Item no. 43154
Black, 45 cm long

PINOFIT® fascia roller Wave

- Hardness: medium
- Extra long 45 cm, ø approx. 12 cm for more exercise options
- Multi-functional training equipment for the most effective exercise
- Specially structured surface with wave profile
- For warming up or cooling down after sport
- Enables multiple exercises with different degrees of difficulty for beginner, advanced and professional levels
- Especially robust inner core for more stability
- Coating of high-quality polyurethane foam for maximum anti-slip properties and shape retention
- Including exercise poster
- With online training video

All colours freely combinable



Item no. 43153
Red, 45 cm long



Item no. 43134
Yellow, 45 cm long



Item no. 43156
Pink, 45 cm long



Item no. 43132
Ice blue, 45 cm long



Item no. 43133
Purple, 45 cm long



Item no. 43199
Deep blue, 45 cm long



Item no. 43196
Grey, 45 cm long



Item no. 43157
Azure, 45 cm long



Item no. 43155
Lime, 45 cm long

PINOFIT® fascia roller Wave Pro

- Hardness: medium
- Approx. 36.5 cm long, ø approx. 13 cm
- The surface with innovative wave profile is designed to intensively stimulate the muscle fibres during use, thus improving blood circulation
- The ergonomically shaped indentation is designed to provide additional relief for the spine
- The multifunctional training tool for a highly effective workout in the practice and at home
- For warming up or cooling down after the workout
- Enables multiple exercises with different degrees of difficulty for beginner, advanced and professional levels
- Especially robust inner core for more stability
- Coating of high-quality polyurethane foam for maximum anti-slip properties and shape retention
- Including exercise poster
- With online training videos

All colours freely combinable



Ergonomically shaped indentation to provide ideal relief for the spine (and other parts of the body) during the workout.



Item no. 43195
Azure, 36.5 cm long




Item no. 43193
Lime, 36.5 cm long



Item no. 43183
Black, 36.5 cm long

FASCIA TRAINING

Diverse tools for various requirements in therapy and training.



PINOFIT® fascia roller Wave Mini

- Hardness: medium
- Optimised for the treatment of feet, hands and arms
- Compact training tool with specially structured surface with wave profile
- Especially robust inner core for more stability
- Coating of high-quality polyurethane foam for maximum anti-slip and shape retention
- Convenient for travelling and suitable for your sports bag

Lime Item no. 43116 Approx. 15 cm long	Black Item no. 43118 Approx. 15 cm long
Azure Item no. 43117 Approx. 15 cm long	Pink Item no. 43119 Approx. 15 cm long

All colours freely combinable



PINOFIT® fascia roller Stick

- Hardness: hard
- Length without handles approx. 28 cm
- The waved surface provides a pleasant massage effect
- Perfect for loosening fascias in arms, legs, core muscles and buttocks
- Special wave-shaped surface
- For gentle pressure control

Lime Item no. 43173 Approx. 28 cm, ø 8 cm	Azure Item no. 43172 Approx. 28 cm, ø 8 cm	Black Item no. 43174 Approx. 28 cm, ø 8 cm
--	---	---

All colours freely combinable




PINOFIT® half ball wave

- Hardness: medium
- Especially for coordination and sensomotoric exercises
- The waved surface provides a pleasant massage effect
- Usable on both sides
- With needle valve for individual pressure regulation

Lime Item no. 43178 Approx. ø 16 cm x 8.5 cm	Azure Item no. 43179 Approx. ø 16 cm x 8.5 cm	Black Item no. 43177 Approx. ø 16 cm x 8.5 cm
---	--	--

All colours freely combinable



PINOFIT® fascia Duo-Ball Wave

- Hardness: hard
- Specially structured surface with wave profile
- Ideal for fascia self-massage and targeted activation of trigger points
- For warming up or cooling down after sport
- Enables multiple exercises with different degrees of difficulty for beginner, advanced and professional levels
- Coating made of high-quality polyurethane foam for maximum anti-slip properties and shape retention
- Practical for on the go

Lime Item no. 43188 Approx. 15 cm x 8 cm	Azure Item no. 43189 Approx. 15 cm x 8 cm	Black Item no. 43187 Approx. 15 cm x 8 cm
---	--	--

All colours freely combinable



PINOFIT® fascia ball Wave

- Hardness: hard
- Convenient for travelling and suitable for your sports bag
- specially structured surface with wave profile
- Ideal for fascia self-massage and targeted activation of trigger points
- For warming up or cooling down after the workout
- Enables multiple exercises with different degrees of difficulty for beginner, advanced and professional levels
- Coating made of high-quality polyurethane foam for maximum anti-slip properties and shape retention

Lime Item no. 43181 Approx. ø 9 cm	Azure Item no. 43182 Approx. ø 9 cm	Black Item no. 43180 Approx. ø 9 cm
---	--	--

All colours freely combinable



HARD FOAM FASCIA ROLLERS

For the highest expectations in therapy, fitness & gym.



45 cm

PINOFIT® fascia roller Foam

- Hardness: hard
- Extra long 45 cm, ø approx. 15 cm and extra light
- Multifunctional training device for effective workout
- Slightly structured surface
- From especially robust polypropylene for the highest level of shape retention
- For warming up or cooling down after the workout
- Permits various exercises with different degrees of difficulty for beginner, advanced and professional levels
- Including exercise poster
- With online training videos



Black
Item no. 43160 extra long 45 cm, ø approx. 15 cm

15 cm

PINOFIT® fascia roller Mini

- Hardness: hard
- 15 cm, ø approx. 5 cm
- Made from especially robust polypropylene
- Ideal for application on small areas, for example soles of feet



Black
Item no. 43166 15 cm, ø approx. 5 cm

42 cm

Ergonomically shaped indentation to provide ideal relief for the spine and other parts of the body during the workout

PINOFIT® fascia roller Pro

- Hardness: medium
- 42 cm, ø approx. 15 cm
- Made from especially robust polypropylene
- Wave profile for pleasant massage effects
- Including exercise poster
- With online training videos



Black
Item no. 43165 42 cm, ø approx. 15 cm

15 cm

PINOFIT® Duoball

- Hardness: hard
- Approx. 15 cm length x 8.5 cm height
- Compact apparatus for the massage of fascias
- Made from very robust polypropylene



Black
Item no. 44295 Approx. 15 cm length x 8.5 cm height

SPIKY MASSAGE BALLS

Perfection all round: price and delivery are on the ball.



16 cm



PINOFIT® fascia spiky roller Peanut


- Compact apparatus for the massage of fascias
- Ideal, for example, for gentle neck massages
- With needle valve for individual pressure regulation

Lime
Item no. 44288
Approx. 16 cm long

Azure
Item no. 44289
Approx. 16 cm long

Dark grey
Item no. 44287
Approx. 16 cm long

All colours freely combinable



PINOFIT® spiky massage ball

- Spiked ball for massage
- Ideal for grip and strengthening exercises
- With needle valve for individual pressure regulation

Yellow semitransparent
Item no. 44833 ø approx. 7 cm

Red semitransparent
Item no. 44208 ø approx. 8 cm

Lime semitransparent
Item no. 44217 ø approx. 9 cm

Azure semitransparent
Item no. 44226 ø approx. 10 cm



TRIGGER & MASSAGE TOOLS

For intensive sensorimotor applications.




Three different degrees of hardness

PINOFIT® spiky massage ball set

- For self-massage and use on trigger points as well as coordination and sensorimotor exercises
- Three different degrees of hardness (soft = pink, medium = lime, hard = azure)
- The ball is pushed back and forth between a hard surface and the body part which needs to be treated

Item no. 44286 3 spiky balls each with ø approx. 6.5 cm



Three different degrees of hardness

PINOFIT® trigger ball set

- For self-massage and use on trigger points
- Three different hardness levels (soft = pink, medium = lime, hard = azure)
- The ball is pushed back and forth between a hard surface and the body part which needs to be treated

Item no. 44285 3 balls each with ø approx. 6 cm



PINOFIT® Balance spiky ball

- Especially for coordination and sensorimotoric exercises
- With needle valve for individual pressure regulation

Lime

Item no. 40521 ø approx. 16 cm

Dark grey

Item no. 40621 ø approx. 16 cm

All colours freely combinable

PINOFIT® fascia spiky massage ball

- Compact apparatus for the sensorimotoric massage of fascias
- With needle valve for individual pressure regulation

Lime

Item no. 44290 Approx. 15 cm long



PINOFIT® Pilates ball

- For strengthening and relaxation exercises
- Useful for Pilates and yoga workouts
- Ideal also for rehab and back strengthening exercise groups
- With plug for a fast and simple change of the ball's volume

Dark grey

Item no. 44225
Approx. ø 18 cm

Lime

Item no. 44228
Approx. ø 18 cm

Red

Item no. 44229 Approx.
ø 22 cm

All colours freely combinable



PINOFIT® exercise ball Antiburst

- Highly elastic special plastic
- Very durable, elastic and safe
- Anti-burst material makes sure that air will escape slowly in case of damage
- Ideal for sports and gymnastics

All colours freely combinable



Item no. 44227
Lime, approx. ø 65 cm



Item no. 44247
Azure, approx. ø 65 cm



Item no. 44237
Dark grey, approx. ø 65 cm



Supplied without products.

Wall mount fixture

- For three balls
- Supplied without fixtures or balls

Item no. 48193 Silver, powder-coated sheet steel, 198 x 45 x 20 cm

Stacking aid

Item no. 45009 1 item = 3 rings (supplied without balls).

Ball base (not shown)

Item no. 45020



Supplied without products.



PINOFIT® bouncy animals

- For active, creative and fun play and recreation
- Perfect toy for the waiting area
- Makes learning balance and motor skills fun
- Suitable for indoors and outdoors
- Made from sturdy flexible plastic
- Easy to clean
- Latex-free
- Can be inflated with standard air pump (not included)
- Sitting height: approx. 25 cm, may vary depending on air volume
- Age: approx. 1.5 years and above
- Maximum load approx. 40 kg

Bouncy animal - dragon

Item no. 38440 L x W x H: approx. 60 x 25 x 51 cm, sitting height: approx. 25 cm

Bouncy animal - giraffe

Item no. 38435 L x W x H: approx. 56 x 29 x 47.5 cm, sitting height: approx. 28 cm

All colours freely combinable

PINOFIT® shoulder stretcher

- For target training and for mobilisation of the shoulder, arm and back muscles
- Easy to use and effective – simply fasten between door and door frame
- The smooth-running pulley ensures light and flowing movement
- for exercises performed while standing or sitting
- Also perfect for home use
- Nylon rope in colour black with two handles



Item no. 43525 Approx. 185 cm long

Traction belt

- Sturdy aircraft webbing with double stitching
- Sturdy metal belt buckle (folding quick release)
- Versatile applications
- Can be individually adjusted



Grey
Item no. 42892
Approx. L x W: 220 x 5 cm



Red
Item no. 42891
Approx. L x W: 220 x 5 cm



PINOFIT® skipping rope

- For warm-up exercise and condition training
- Trains speed, coordination, jumping power and endurance
- Perfect for exercise at the practice or at home
- Steel rope with soft coating with two anti-slip and ergonomically shaped handles

Item no. 44658 Length (without handles), adjustable from approx. 268 cm to 296 cm

PROFESSIONAL YOGA MATS

Especially non-slip and dimensionally stable thanks to PU underside.



With positioning print.

PINOFIT® yoga mat with positioning print

- Innovative PU surface with positioning print makes it easier to assume the starting position
- Top and bottom in two contrasting colours
- Non-slip
- Roll up to save space and for easy transport
- Ultra light, weight approx. 1.2 kg

Lime/grey

Item no. 42858 Approx. 182 x 61 x 0.4 cm

Fuchsia/grey

Item no. 42859 Approx. 182 x 61 x 0.4 cm

All colours freely combinable

PINOFIT® yoga mat

- Mat with non-slip under and upper for a firm grip
- Non-elastic, for a firm grip
- Top and bottom in two contrasting colours
- Roll up to save space and for easy transport
- Ultra light, weight approx. 0.6 kg

Lime/grey

Item no. 42851 Approx. 183 x 61 x 0.4 cm

Warm grey/grey

Item no. 42868 Approx. 183 x 61 x 0.4 cm

All colours freely combinable

PINOFIT® reflex zone mat

- Perfect for stimulating and massaging the soles of the feet and back area
- Rollable & space-saving
- Lightweight, durable EVA material

Silver-grey

Item no. 42838 Approx. 120 x 35 x 1.8 cm



NON-SLIP EXERCISE MATS

Top quality at fair prices.



Dark grey/black

Item no. 42842 Approx. 186 x 61 x 1.0 cm

Lime/black

Item no. 42850 Approx. 186 x 61 x 1.0 cm

Ice blue/black

Item no. 42841 Approx. 186 x 61 x 1.0 cm

All colours freely combinable

PINOFIT® exercise mat

- Mat with non-slip bottom and top for a firm grip
- High absorption
- Top and bottom in two contrasting colours
- Roll up to save space and for easy transport
- Ideal for fitness and gymnastic exercises and much more
- Weight approx. 1.5 kg

PINOFIT® mat trolley

- For 10–20 mats, 60 cm wide
- Powder-coated sheet steel
- Extremely robust and durable
- Low rolling resistance wheels for mobile use
- Supplied unassembled and without products

Item no. 48281 W x D x H: 132 x 70 x 107 cm



PINO wall mount fixture


- For 5–10 mats, supplied without fixtures or mats
- Silver powder-coated sheet steel

Item no. 48194 W x D x H: 135 x 25 x 20 cm



BALANCE PADS

The classic for effective balance training.



PINOFIT® balance pad

- Multi-purpose workout tool
- Used for balance and coordination exercises
- Workout tool for more stability and muscle strength
- Including exercise poster

Lime

Item no. 44545 Approx. 50 x 39 x 6 cm

Azure

Item no. 44555 Approx. 50 x 39 x 6 cm

Dark grey

Item no. 44560 Approx. 50 x 39 x 6 cm

All colours freely combinable



MEDICAL PRODUCT CE



AIREX® balance pad

- Because of its destabilising properties, can be used for training motor skills in rehabilitation, prevention and competitive sports
- For combined training to improve conditional and coordinative skills
- Closed-cell, super-soft special foam
- Does not absorb moisture or dirt

Item no. 44502 Approx. 50 x 41 x 6 cm



AIREX® balance pad Elite

- Same effect as balance pad but more slip-resistant due to ribbed surface on both sides, and therefore more durable
- For use as a therapy and workout device to restore walking or foot function

Blue

Item no. 44503 Approx. 50 x 41 x 6 cm

Lava


Item no. 44506 Approx. 50 x 41 x 6 cm

THERAPY MATS

The Swiss original for 50 years.

- Perfect for remedial and physical therapy, physiotherapy and rehabilitation, recreational and apparatus gymnastics, yoga, sports for the disabled and much more
- Closed cell, therefore does not absorb dirt or moisture
- Does not cause friction-induced burns, is lightweight, buoyant, resistant to aging, insulates body temperature and has optimal damping properties

MEDICAL PRODUCT CE



AIREX® Coronella 200

Approx. 2.5 kg weight.

Slate

Item no. 45935 Approx. 200 x 60 x 1.5 cm

Platinum

Item no. 45940 Approx. 200 x 60 x 1.5 cm



AIREX® Coronella 185

Approx. 2.3 kg weight.

Red

Item no. 45901 Approx. 185 x 60 x 1.5 cm

Slate

Item no. 45919 Approx. 185 x 60 x 1.5 cm

Blue

Item no. 45917 Approx. 185 x 60 x 1.5 cm

AIREX® Hercules 200 (not shown)

Approx. 7.2 kg weight.

Blue

Item no. 45956 Approx. 200 x 100 x 2.5 cm

Green

Item no. 45947 Approx. 200 x 100 x 2.5 cm



AIREX® Corona 185

Approx. 3.9 kg weight.

Red

Item no. 45929 Approx. 185 x 100 x 1.5 cm

Slate

Item no. 45937 Approx. 185 x 100 x 1.5 cm

Green

Item no. 45938 Approx. 185 x 100 x 1.5 cm

Blue

Item no. 45939 Approx. 185 x 100 x 1.5 cm

AIREX® Fitline 180

Approx. 1.8 kg weight.

Kiwi

Item no. 45908 Approx. 180 x 60 x 1.0 cm

Slate

Item no. 45909 Approx. 180 x 60 x 1.0 cm

Water blue

Item no. 45905 Approx. 180 x 60 x 1.0 cm

AIREX® Atlas (not shown)


Approx. 5.9 kg weight.

Red

Item no. 45885 Approx. 200 x 125 x 1.5 cm

Green

Item no. 45894 Approx. 200 x 125 x 1.5 cm



AIREX® balance beam

- For balance training, for agility, balance and posture stability training
- For coordination and response training and for "Balanced Motion" endurance training


Item no. 45930 Approx. 162 x 24 x 6 cm



SoftX® Pilates roller

- Optimal training device for the trend sport Pilates
- For efficient training at home or group workouts at the gym
- Slip-resistant, with good grip
- Skin-friendly surface

Item no. 45204 Approx. Ø 14.5 cm, length approx. 90 cm




AIR MAT

A revolution in greater dynamics.



Foot pump
For air mat & air pad
Item no. 44340



Inner drop-stitch technology for maximum stability.



Different dynamics possible.



With practical carrying handle.

PINOFIT® air mat

- Air-filled mat
- Vary the training intensity by adding or removing air
- Maximum stability thanks to two special plastic layers interwoven with nylon threads (drop-stitch technology)
- Disinfectable, stable and non-slip surface
- Printed markings for correct posture
- Delivery incl. adapter for quick inflation with standard foot pump (not included)
- Can be rolled up to save space
- Max. load 200 kg

Item no. 44330 Grey, approx. 185 cm x 85 cm x 10 cm

AIR PAD

The compact complement to the air mat.



PINOFIT® air pad

- The more dynamic balance pad
- Space-saving version with the same properties as the PINOFIT® air mat
- Perfect for targeted training sessions and group workouts
- Max. load 200 kg

Item no. 44335 Grey, approx. 50 cm x 50 cm x 10 cm



Practical safety screw valve.




Control the dynamics simply by removing air.

PINOFIT® sports mat active physio mat

- Facilitates communication in individual and group training
- Suitable for documentation of movement dimensions
- Punched division into fields with centre circle
- Non-slip underside with wave structure

Item no. 42843 Grey, approx. 150 x 100 cm



HAND THERAPY

For varied & individual training.

PINOFIT® hand trainer

- Continuously adjustable resistance
- Small, light & compact: ideal for on the go
- Ergonomically shaped
- Anti-slip handles for stable grip
- Stainless steel spring
- The strength is based on the amount of resistance when the trainer is squeezed


Yellow

Item no. 41160 10–40 kg, L x W x D: 15 x 10.5 x 2.3 cm

Grey

Item no. 41165 10–60 kg, L x W x D: 15 x 10.5 x 2.3 cm

All colours freely combinable



Handmaster plus

- Exercises the stretching and contraction of the relevant hand muscles
- The training is intended to support the improvement of the hand and arm muscle coordination

Purple

Item no. 41110 Soft

Red

Item no. 41113 Medium

Orange

Item no. 41115 Strong

All colours freely combinable



PINOFIT® fitness clay

- High-quality fitness clay
- for exercising finger and hand muscles

Yellow

Item no. 45320 Extra soft, 85 g

Pink

Item no. 45321 Soft, 85 g

Lime

Item no. 45322 Medium-soft, 85 g

Blue

Item no. 45323 Firm, 85 g

All colours freely combinable



New



PINOFIT® training ball with positioning splint

- Stabilises the wrist thanks to removable splint
- Adjustable wrist cuff & fingers to customise size
- With hook-and-loop fastener

Grey

Item no. 41144 Extra soft, 80 g, approx. 9 cm

New



PINOFIT® hand trainer

- Ergonomic hand trainer for finger and hand exercises

Grey

Item no. 41142 Approx. 11 x 10 cm



New



PINOFIT® hand trainer Flex

- Resistance variably adjustable for each finger
- With spring-loaded keys for training individual fingers
- Sliding button for special training of the thumb
- Approx. 1.1–2.7 kg, individually adjustable resistance

Grey

Item no. 41140 Resistance can be adjusted individually from 1.1 to 2.7 kg

New



Set of 4

PINOFIT® Squeeze Egg set of 4

- Different resistance levels in different colours
- Ergonomically formed to fit shape of hand

Item no. 41175 Yellow = light resistance, 15 kg, red = medium resistance, 20 kg, green = strong resistance, 25 kg, blue = extra strong resistance, 30 kg

FOOT TRAINING

Intensive stretching & relaxing.



Set of 2

Foot-stretcher

- Floor protectors prevent slipping and ensure a gentle foot movement
- Rubber stopper guarantees stability
- With heightened edges for a secure hold
- Made from robust polypropylene

Black

Item no. 43527 2 pieces, approx. 31.5 cm x 8.5 cm



PINOFIT® Separatoes

- Stretching of foot structures
- Use for sport activities such as yoga and Pilates
- Supports the relaxation of feet and legs after a long day
- Made from high quality flexible silicone
- Flexible and soft
- Easy to clean
- Latex-free

Set of 2

Light blue

Item no. 42860 1 pack small, 1 pair

Set of 2

Deep blue

Item no. 42861 1 pack large, 1 pair



DUMBBELL WITH SILICONE COATING

Particularly good grip & easy to disinfect.



New



Set of 2

PINOFIT® weight cuff

- With soft touch silicone sheath
- Ergonomic design
- Easy to clean with alcohol-free disinfection

Dark grey

Item no. 44751 0.5 kg

Dark grey

Item no. 44752 1 kg

New



Set of 2

PINOFIT® ring dumbbell with handle

- With soft touch silicone sheath
- Ergonomic design
- Easy to clean & disinfect

Dark grey

Item no. 44757 0.5 kg, 2 pieces

Dark grey

Item no. 44758 1 kg, 2 pieces

Dark grey

Item no. 44759 2 kg, 2 pieces

New



Set of 2

PINOFIT® barbell

- With soft touch silicone sheath
- Ergonomic design
- Easy to clean & disinfect

Dark grey

Item no. 44754 0.5 kg, 2 pieces

Dark grey

Item no. 44755 1 kg, 2 pieces

Dark grey

Item no. 44756 2 kg, 2 pieces

ROCKER BOARDS

For multi-functional training.



PINOFIT® rocking board

- Ideal for balance and posture exercises in sitting, lying or kneeling positions
- For balancing and rocking
- Stable construction made from birch plywood, natural
- With anthracite, anti-slip top
- Maximum load 100 kg
- Rubber lining for floor protection



Feet with non-slip rubber.



Item no. 41265 Approx. 50 x 50 cm



Item no. 41260 Approx. 90 x 50 cm

BALANCE BOARDS MADE FROM WOOD

The sustainable training alternative.



With medium movement dynamics

PINOFIT® balance board

- For training to improve the stability of the leg axis, balance and coordination
- Large platform with non-slip coating for hip-width exercises
- Also for balance exercises sitting or lying down or as arm support
- Birch multiplex, natural
- One damper
- Maximum load 100 kg
- Non-slip thanks to rubber studs

Item no. 44541 Approx. ø 49 cm

With low movement dynamics

PINOFIT® balance board Rehab

- For training to improve the stability of the leg axis, balance and coordination
- Large platform with non-slip coating for hip-width exercises
- Also for balance exercises sitting or lying down or as arm support
- Birch multiplex, natural
- One damper
- Maximum load 100 kg
- Non-slip thanks to rubber studs

Item no. 44542 Approx. ø 49 cm

HIGH-QUALITY BALANCE CUSHION

For a strong back and better balance.



Incl. ball pump.



With needle valve for individual pressure adjustment

PINOFIT® balance cushion ø 33 cm

- Functional and versatile exercise tool for balance and strength-building exercises
- For back-friendly seating, for example on your office chair
- One side smooth and the other side with soft nubs
- The tissue of the area that is being exercised can be additionally stimulated by using the nub-studded side
- Including ball pump and exercise poster

Azure

Item no. 43136 Approx. ø 33 cm

Deep blue

Item no. 43144 Approx. ø 33 cm

Lime

Item no. 43137 Approx. ø 33 cm

Dark grey

Item no. 43142 Approx. ø 33 cm

All colours freely combinable



Therapy balance pad

- Therapy balance pad for balance and coordination exercises
- Robust, Easily disinfected PVC quality

Item no. 44530 Soft green, approx. ø 39 cm

TOGU® Jumper®

- Ideal pressure distribution thanks to proven Dynair® technology
- Extremely dynamic rebound (trampoline effect)
- More stability for ankle joints
- Non-slip on all surfaces
- Also designed for lateral jumping movements
- Better grip surface
- Usable on both sides
- Vibrations transferred to the body
- Delivery includes exercise poster

Item no. 44300 Red, approx. ø 52 cm

HIGH-QUALITY BALANCE TRAINER

Dynamic exercise support.




PINOFIT® balance trainer

- Half ball from highly elastic special plastic
- Multi-functional training tool for balance and coordination exercises
- The dynamic trampoline effect transfers vibrations to the body
- Usable on both sides
- Pressure can be individually adjusted
- Maximum load 120 kg



Incl. ball pump.



PINOFIT® balance trainer

Dark grey/black

Item no. 43127 Approx. ø 60 cm, max. 22 cm tall

Dark grey/lime

Item no. 43128 Approx. ø 60 cm, max. 22 cm tall

All colours freely combinable



PINOFIT® balance trainer Compact

Dark grey/black

Item no. 43122 Approx. ø 47 cm, max. 16 cm tall

Dark grey/lime

Item no. 43123 Approx. ø 47 cm, max. 16 cm tall

All colours freely combinable

SPECIAL PRODUCTS

For varied training.



PINOFIT® Physio Chair XL

- Versatile for exercises in sitting and lying position
- Copes with all demands in your daily practice for therapy, rehab and simple fitness training
- As a balance stool for dynamic sitting – the gentle rocking of the Physio Chair promotes active sitting correction and supports an upright sitting position
- As a balance trainer for improving balance and coordination
- As a fascia roller for warming up or cooling down after sport
- As a yoga wheel for stretching torso and back muscles

Lime

Item no. 42867 Height approx. 52 cm, ø approx. 30 cm

Dark grey

Item no. 42869 Height approx. 52 cm, ø approx. 30 cm

All colours freely combinable

DUMBBELLS AT TOP PRICES

Very easy to hold.



PINOFIT® kettlebell PVC Soft & Safe

- For a versatile, functional whole-body workout
- Helps to strengthen the deep-seated muscle groups
- Lower risk of injuries
- Ergonomic, dimensionally stable handle
- Elastic, soft underside protects the floor
- Made from PVC plastic
- Filled with sand
- Available in different weight variants

2 kg

Item no. 42952 Dark grey

4 kg

Item no. 42954 Dark grey

6 kg

Item no. 42956 Dark grey



Reebok® Step

- Multi-purpose workout tool
- With anti-slip rubber step area and rubber feet for stable footing
- Height can be adjusted to 15, 20 or 25 cm with the help of mono blocks
- Maximum load 110 kg

Item no. 42837

Approx. 90 x 35 cm, weight approx. 7.5 kg



PINOFIT® Physio Wheel

- Ideal for stretching and balance exercises
- Supports, facilitates and varies workouts
- Non-slip, disinfectable surface
- Made from special dimensionally stable plastic
- Max. load 150 kg

Item no. 42890

Black approx. 39 cm x 20 cm x 13 cm



PINOFIT® hula hoop

- The fun way to train balance and coordination
- Perfect for back exercises
- Eight elements can be connected to each other easily thanks to the click system
- Can be taken apart quickly to save space
- Incl. storage bag

Grey/lime

Item no. 44666 Approx. ø 96 cm

Grey/black

Item no. 44667 Approx. ø 96 cm

All colours freely combinable



Set of 2

PINOFIT® dumbbells

- Very easy to hold
- Easy to clean and disinfect
- Soft PVC coating protects floors and mats
- Sweat- and scratch-resistant
- 0.5 to 5 kg available

0.5 kg

Item no. 44770
Dark grey, 2 pieces

1 kg

Item no. 44772
Dark grey, 2 pieces

1.5 kg

Item no. 4477
Dark grey, 2 pieces



2 kg

Item no. 44774
Dark grey, 2 pieces

3 kg

Item no. 44775
Dark grey, 2 pieces

4 kg

Item no. 44776
Dark grey, 2 pieces

5 kg

Item no. 44777
Dark grey, 2 pieces




PHYSIO BOOST MINI

The massager in pocket format.



Discover our compact mini massager.



The Physio Boost is easy to transport in its practical bag, always ready to hand and quick to use. I like to use it as a support for wellness and well-being massages – and not only because it makes massaging easier on my hands. The main benefit is that the vibrations can have a positive effect on the tone of the muscles.

Carsten Boss, Senior Therapist



PINOFIT® Physio Boost Mini

- Compact mini massager with non-slip wristband
- Lightweight, ergonomic and cordless design
- Powerful lithium-ion battery (2,500 mAh)
- Convenient charging via USB
- Particularly quiet (approx. 50 dB)
- 4 intensity levels with 900 to 3,200 massage strokes per minute
- 4 interchangeable massage attachments for different applications:
 - Small, round attachment for full-body massages and small muscle groups
 - Flat attachment for full-body massages and stimulating massages
 - Bullet attachment with removable silicone damper for muscle and trigger point massages
 - Double-head attachment with removable silicone dampers for muscle massages
- Battery runtime: approx. 7 hours
- Max. charging time: approx. 8 hours via USB
- Approx. 142 x 90 x 44 mm, approx. 512 g
- Please observe the safety instructions in the operating information to avoid personal injury and damage to property.

Silver
Item no. 44696


Petrol
Item no. 44694

Kiwi
Item no. 44697

Black
Item no. 44693

All colours freely combinable

NEW
New colours



Four interchangeable massage attachments: round, double-head, flat, bullet.



With practical storage bag.

SILICONE TRIGGER TOOLS

For spot and area application.



PINOFIT® Physio Tool

- Different tool sides for triggering, mobilising, massaging and stroking out
- Makes work easier for the therapist's hands
- Made from high-quality silicone
- Flexible and soft
- Multi-purpose
- Durable and shatter-proof
- Easy to clean
- Latex-free

Item no. 35814 Dark grey, L x W x H, approx. 13 x 8 x 4 cm

PINOFIT® trigger point & massage tool

- Ideal for massages and for use on trigger points
- For the supporting treatment of fascial tensions
- Also useful for the self-massage at home
- Made from high-quality and durable silicone

Coral Item no. 39317 Soft Lime Item no. 39316 Medium Dark grey Item no. 39315 Hard

All colours freely combinable

PINOFIT® trigger point & massage tool 3er Set

Item no. 39318 Set of three tools: 1x coral (soft), 1x lime (medium), 1x dark grey (hard)



ALUMINIUM FASCIA TOOLS

For instrument-supported mobilisation.



PINOFIT® fascia tools for modern mobilisation supported with instruments (IASTM: instrument-assisted soft tissue mobilisation) facilitate the intensive treatment of dysfunctions caused by fascial tension and adhesion. Therapeutic knowledge and manual techniques are preconditions for a successful application. Following intensive research into functions, PINOFIT® fascia tools serve to improve the

effects of several hands-on-techniques when treating superficial and deeper fascias and retinacula.

All tools have a bevelled edge. This increases the contact surface. Changing the angle of application thus creates a change in pressure intensity.

✓ MEDICAL PRODUCT CE



Fascia tool large for large body regions.



Fascia tool small for small body regions.



Fascia tool large

- The concave or straight sides of the tool are used specifically for large body and fascia areas
- It is possible to stretch fascia with the spoon-shaped end of the tool while providing a secure hold for the thumb

Item no. 48212 Hand-made from powder-coated aluminium, length approx. 27 cm



Fascia tool small

- Concave tool side for the special treatment of small body and fascia areas
- Furthermore, the paravertebral treatment and the fixation of vertebrae can be carried out with the wave-shaped tool side with a minimum of effort

Item no. 48211 Hand-made from powder-coated aluminium, length approx. 21 cm



Set of 2

Fascia therapy set

- 1 fascia tool large
- 1 fascia tool small
- With online training videos

Item no. 43164

MASSAGE CUPS

Effective support in treatment.



Set of 2

Massage cups Soft Touch

- for wellness massages
- Can easily be moved on the skin
- Made from high-quality silicone
- Flexible and soft
- Multi-purpose
- Durable and shatter-proof
- Easy to clean
- Latex-free
- Extra-soft surface

Ice blue

Item no. 35485 2 pieces, ø 5 cm



Soft & flexible

Massage cups

- for wellness massages
- Can easily be moved on the skin
- Made from high-quality silicone
- Flexible and soft
- Multi-purpose
- Durable and shatter-proof
- Easy to clean
- Latex-free

Set of 2

Item no. 35470 Small, 2 pieces, ø 3.5 cm

Item no. 35471 Medium, 2 pieces, ø 5 cm

Item no. 35472 Large, 1 piece, ø 7 cm



Set of 2

Massage cup set

Item no. 35474 2 pieces small, approx. ø 3.5 cm, 2 pieces medium, approx. ø 5 cm 1 piece large, approx. ø 7 cm

FASCIA CUPS

For professional massage and gentle cupping treatments.



Fascia cups

- Ideal for spot massage techniques
- Simple creation of a vacuum, without any further aids
- Made from high-quality silicone
- Extra-strong suction effect
- Durable and shatter-proof
- Easy to clean
- Latex-free



Item no. 35480 Small, 2 pieces approx. ø 2.5 cm

Set of 2

Item no. 35481 Medium, 2 pieces approx. ø 5.0 cm

Set of 2

Item no. 35482 Large, 1 piece approx. ø 6.5 cm



Fascia cup set

Item no. 35483 2 pieces small, approx. ø 2.5 cm, 2 pieces medium, approx. ø 5.0 cm, 1 piece large, approx. ø 6.5 cm



Set of 6

Trigger cups

- Ideal for spot massage techniques
- For targeted cupping on trigger points
- Especially wide and flat shape for easy application
- Simple creation of a vacuum, without any further aids
- Made from high-quality silicone with inlaid tourmaline quartz pearls
- Durable and shatter-proof
- Easy to clean
- Latex-free

Item no. 35490 Pack with 6 pieces, each approx. ø 6 cm



Set of 6



Trigger cups Pro

- Ideal for spot massage techniques
- For trigger point and cupping application
- Particularly easy to grip and handle
- Practical helpers that make work easier for your hands
- Especially wide and flat shape for easy application
- Simple creation of a vacuum, without any further aids
- Made from high-quality silicone
- Extra-strong suction effect
- An alternative to traditional glass suction cups
- Durable and shatter-proof
- Easy to clean
- Latex-free

Item no. 35475 Pack with 6 pieces, each approx. ø 4 cm



Trigger stick

- For targeted and optimally dosed applications on trigger points
- Made from solid wood
- For support of massage techniques
- Two treatment points with different profiles
- Firm grip
- Fits into any trouser pocket
- Also suitable for self-treatment and for use by patients

Item no. 36507

NECK PILLOWS

For maximum sleeping comfort.

Perfect if you sleep on your stomach

DREAM ME! pillow Butterfly

- Perfect heat regulation thanks to a core of viscose elastic memory foam and special gel applications for a pleasant sleep ambience
- Due to the special butterfly shape with gaps on the sides suitable for all shoulder positions and in particular for people who sleep on their stomach
- For back, stomach and side position
- The ergonomic shape is supportive, guarantees a perfect pressure distribution and adapts perfectly to the contours of your head and neck
- Cover is washable at 40°C

Item no. 45266 Approx. 58 x 36 x 7/11 cm



Adjustable height



Viscoelastic high-tech material.



- Sleep like a baby and recharge your batteries! Wake up in the morning with renewed energy, and enjoy a day full of get-up-and-go!
- Hi-tech viscoelastic material developed in space research
- The wave shape is ergonomically supportive and fits perfectly to the contours of your head and neck
- For back and side position
- With additional variable inserts for individual adjustment of the height to your shoulder width
- The soft, removable comfort cover made of velour and the rigid viscoelastic core ensure long-lasting quality
- Velour cover washable at 40°C
- Adjustable height

DREAM ME! pillow

Item no. 45286 Incl. protective cover, approx. 50 x 30 x 11/10/8/7 cm

DREAM ME! pillow XL

Item no. 45276 Incl. protective cover, approx. 60 x 40 x 13/12/11/10 cm

DREAM ME! pillow replacement cover

Item no. 45288 Approx. 50 x 11 x 30 cm



With non-slip underside.



Cushion Special

- Polyurethane core
- The shaped leg parts reduce the pressure on the coccyx and support an upright posture to enable relief of the spine
- Underside with anti-slip nubs
- Cover washable at 30°C (delicate cycle)
- Breathable, with zip
- Ideal for office chairs or for travel

Item no. 44663 Approx. 42 x 40 x 9/2cm



With special gel chambers.



Memory foam cushion

- Perfect heat regulation and pleasant seating climate thanks to a core of viscose elastic memory foam and special gel chambers
- The ergonomic shape adapts perfectly to any form of seating
- The U-shape with flexible, perfectly shaped leg parts reduces the pressure on the coccyx and supports an upright posture to enable relief of the spine
- Underside with anti-slip nubs
- Cover washable at 30°C (delicate cycle)
- Breathable, with zip
- Ideal for office chairs or for travel

Item no. 44647 Approx. 46 x 37 x 8 cm

KINESIOLOGY TAPE



PINOTAPE
THE THERAPIST'S CHOICE
SINCE 2006



PINOTAPE THE TAPE OF THERAPISTS




Reliable top quality
for professional applications.

The right tape for every patient!



TOP
SELLER



	PINOTAPE Sport	PINOTAPE Therapy	PINOTAPE Sensitive
Material	 Rayon (viscose) <ul style="list-style-type: none"> • Quick-drying • Water- and sweat-resistant • Very comfortable to wear • High-performance tape 	 Cotton <ul style="list-style-type: none"> • Natural fibre • Water- and sweat-resistant • Very comfortable to wear 	 Rayon or cotton
Adhesive	Very good skin compatibility	Very good skin compatibility	Without rosin especially skin friendly
Adhesiveness	●●●●●	●●●●○	●●●○○ Extra gentle
Approx. wearing time* up to	5-7 days	5-7 days	1-3 days 2-12 hours for beauty tape
Ideal for	<ul style="list-style-type: none"> • All taping applications in physiotherapy In particular: <ul style="list-style-type: none"> • with long application duration • for intensive athletic strain • For extreme exposure to stress 	<ul style="list-style-type: none"> • All taping applications in physiotherapy • Also ideal as an everyday tape 	<ul style="list-style-type: none"> • All taping applications in physiotherapy In particular: <ul style="list-style-type: none"> • For sensitive skin • For tape applications on sensitive parts of the body (e.g. crook of the arm, face) • Also suitable for people with rosin allergy
Available sizes	5 cm x 5 m the standard 5 cm x 31 m for frequent users 7.5 cm x 5 m for lymph tapes 10 cm x 5 m for lymph tapes on large areas of the body 2.5 cm x 5 m for tape applications on the hand	5 cm x 5 m	5 cm x 5 m 2.5 cm x 5 m for taping on the face (cosmetic taping) Beauty taping

* Average adhesion of the tape in majority of tape applications. The wearing time is significantly influenced by the part of the body that is taped (varying degrees of mechanical stress), the application technique (stronger tension leads to a shorter wearing time) and other factors (e.g. condition of the skin, sporting activity).

ORIGINAL PINOTAPE SPORT

Extra strong adhesion and quickly drying rayon.



PINOTAPE
THE THERAPIST'S CHOICE
SINCE 2006

✓ MEDICAL PRODUCT CE

LABORATORY-TESTED

PINOTAPE® Sport

- Carrier material made of quick-drying rayon (viscose)
- Designed for universal application
- Can withstand all sport activities
- Adhesive applied in wave form with extra strong adhesion for high exposure to stress
- Emulates the properties of human skin
- Suitable even for patients with excessive perspiration
- Very good skin compatibility
- Very long wearing period
- Expansion 140-150%
- Latex-free
- 5 m x 5 cm

You can find our pre-cut strips online at pinoshop.com

All colours and 5-metre varieties freely combinable

Item no. 45022 Blue	Item no. 45121 Ice blue	Item no. 44090 Deep blue	Item no. 45169 Olive	Item no. 45167 Blue diamond	
Item no. 45149 Petrol	Item no. 45151 Mint	Item no. 45026 Lime	Item no. 45091 Red	Item no. 45083 Purple	
Item no. 45117 Coral	Item no. 45021 Orange	Item no. 45092 Yellow	Item no. 45081 Light beige	Item no. 45089 Black	Item no. 45109 Dark grey
Item no. 45086 Silver	Item no. 45098 Green flower	Item no. 45093 Camouflage	Item no. 45126 Anchor	Item no. 45124 Butterfly	Item no. 45128 Dots
Item no. 45153 Savannah	Item no. 45154 Pepita	Item no. 45156 Reptile	Item no. 45157 Jolly roger blue	Item no. 45158 Jolly roger pink	Item no. 45159 Pride



Extra wide for large muscle groups: 7.5 cm & 10 cm

PINOTAPE® Sport XL & XXL

- For more therapy options
- Enables lymph drainage in the case of large-scale blunt injuries
- For shoulder girdle corrections
- Ideal for ligament techniques on large joints

XL light beige

Item no. 45096 5 m x 7.5 cm

XXL light beige

Item no. 45097 5 m x 10 cm



NEW
Especially for
the hands

Extra narrow for the hands: 2 rolls of 2.5 cm x 5 m

PINOTAPE® Sport Hands

- Ideal for taping wrists and fingers
- For stress in everyday life and sports
- Provides extra stability
- Can help to increase mobility in the wrist joint
- Latex-free
- 2 rolls, each 5 m x 2.5 cm

Item no. 45182 Dark grey, 2 rolls, 5 m x 2.5 cm



PINOTAPE® Sport Fascia

- Developed especially for the use in fascial therapies
- Drop-shaped embossing guarantees more flexibility
- Can withstand all sport activities
- 5 m x 5 cm

Item no. 45088 Deep blue



All colours and 5-metre varieties freely combinable

PINOTAPE Sport 31 m roll

- Ideal for clinics & sports practices
- 31 m x 5 cm

All colours freely combinable for 31 metres

Item no. 45164 Petrol	Item no. 45166 Pink	Item no. 45163 Red	Item no. 45137 Blue	Item no. 45136 Lime	Item no. 45138 Light beige	Item no. 45135 Black



ORIGINAL PINOTAPE SENSITIVE

Especially gentle for highly sensitive skin.

✓ MEDICAL PRODUCT CE  LABORATORY-TESTED 

PINOTAPE® Sport Sensitive

- Carrier material made of quick-drying rayon (viscose)
- Extra mild adhesive for very good skin tolerance
- Designed for universal application
- Adhesive applied in wave form for proven bond strength
- Emulates the properties of human skin
- Medium wearing time
- Expansion 140-150%
- Latex-free
- 5 m x 5 cm

New	New	New
Item no. 45171 Light blue	Item no. 45174 Light beige	Item no. 45172 Light black
New	New	New
Item no. 45176 Light lime	Item no. 45173 Light red	Item no. 45178 Shooting star
New		
Item no. 45177 Zoo	Item no. 45129 Plain	



NEW
Now in new colours and designs



PINOTAPE® Therapy Sensitive

- Extremely mild adhesive
- Very good skin compatibility
- Carrier material made from natural cotton
- Reliable adhesion
- Designed for universal application
- Long-term wearability
- Expansion 140-150%
- Latex-free
- 5 m x 5 cm

Item no. 45029 Plain



Sensitive for sensitive skin

⊕ Extremely mild adhesive

All colours and 5-metre varieties freely combinable

ORIGINAL PINOTAPE BEAUTY

Gentle kinesiology taping for facial treatments.



NEW
Especially for the face

PINOTAPE® Sensitive Beauty

- Ideal for taping on the face, neck and décolleté
- Carrier material made of light rayon (viscose)
- Extra mild adhesive for very good skin tolerance
- Adhesive applied in wave form for proven bond strength
- Emulates the properties of human skin
- Expansion 140-150%
- Latex-free
- 2 rolls, each 5 m x 2.5 cm

Item no. 45179 pink, 2 rolls, 5 m x 2.5 cm

All colours and 5-metre varieties freely combinable



FISKARS®

FISKARS® special scissors

- For easy cutting of PINOTAPE® Sport & Therapy
- Cuts exact, clean edges, saving time during application
- Top quality made from stainless steel by Fiskars®

Item no. 33190



PINOTAPE® Display

- Suitable for the presentation of a PINOTAPE of your choice (5 m x 5 cm)
- With free space for individual price communication
- Stable thanks to sturdy construction

Item no. 45175 Set up approx. 23 x 25 x 10 cm for all 5 m x 5 cm PINOTAPES

ORIGINAL PINOTAPE THERAPY


Kinesiology tape for short to medium wearing time.

✓ MEDICAL PRODUCT CE

LABORATORY-TESTED

PINOTAPE® Therapy

- Carrier material made from cotton
- Designed for universal application
- Adhesive applied in wave form for proven bond strength
- Emulates the properties of human skin
- Very good skin compatibility
- Medium wearing time
- Expansion 140-150%
- Latex-free
- 5 m x 5 cm



All colours and 5-metre varieties freely combi-

ORIGINAL PINOTAPE® HEMP

With natural hemp fibres.

✓ MEDICAL PRODUCT CE

LABORATORY-TESTED

PINOTAPE® Hemp

- Cotton with sustainable hemp fibres
- Universal use, high-quality, reliable
- Reliable adhesion
- Very good skin compatibility
- Long-term wearability
- Latex-free
- 5 m x 5 cm

Item no. 45148 Plain



Hemp: A material with a future

- ⊕ Sustainable resource
- ⊕ Quickly renewable
- ⊕ Universal application

All colours and 5-metre varieties freely combi-



PINOTAPE® Prepare Spray

- The professional, highly concentrated spray improves the adhesiveness of the PINOTAPE®
- The skin is perfectly cleansed, refreshed and prepared for the subsequent application

Item no. 31748 50 ml pump spray

PINOTAPE® Remover Spray

- Reliably and easily removes tape bandages and adhesive residue from the skin
- Spray directly onto the tape dressing, allow to act briefly
- Carefully remove tape, wipe off residue


Item no. 31480 75 ml can

All varieties freely combinable




PINOFIT® CROSS TAPE

Perfect for trigger points.



Three colours available



✓ MEDICAL PRODUCT CE

PINOFIT® Cross Tape

- PINOFIT® Cross Tape is applied to the trigger points (also known as pain, reflex and acupuncture points) that were determined manually by the treating therapist
- Breathability and moisture transport are barely impeded at all by the lattice structure of the patch, and the skin continues to be physiologically balanced
- To use on its own or in combination with PINOTAPE®

All colours & sizes freely combinable

Light beige, small

Item no. 49853 1 pack with 20 sheets with 9 patches each (2.7 cm x 2.1 cm)

Light beige, medium

Item no. 49854 1 pack with 20 sheets with 6 patches each (3.5 cm x 3 cm)

Light beige, large

Item no. 49855 1 pack with 20 sheets with 2 patches each (5.2 cm x 4.5 cm)

Red, small

Item no. 49863 1 pack with 20 sheets with 9 patches each (2.7 cm x 2.1 cm)

Red, medium

Item no. 49864 1 pack with 20 sheets with 6 patches each (3.5 cm x 3 cm)

Red, large

Item no. 49865 1 pack with 20 sheets with 2 patches each (5.2 cm x 4.5 cm)

Blue, small

Item no. 49873 1 pack with 20 sheets with 9 patches each (2.7 cm x 2.1 cm)

Blue, medium

Item no. 49874 1 pack with 20 sheets with 6 patches each (3.5 cm x 3 cm)

Blue, large

Item no. 49875 1 pack with 20 sheets with 2 patches each (5.2 cm x 4.5 cm)

CLASSIC SPORT TAPES

Top quality at reasonable prices.

NEW
Improved
quality

✓ MEDICAL PRODUCT CE

PINOFIT® sport tape classic

- For classic tape bandages
- Non-elastic adhesive tape
- High adhesiveness
- Easy to tear lengthwise and across
- Easy to unroll
- 100% cotton
- Water-repellent
- Latex-free

Item no. 45101 White, 3.8 cm x 10 m



Our highlights
for your sports tape

- ⊕ better adhesive properties
- ⊕ Even easier to tear lengthwise and across

Self-adhesive



✓ MEDICAL PRODUCT CE

PINOFIT® adhesive tape

- Support bandage for producing light compression as an initial measure for minor sports injuries
- Elastic, self-adhesive
- 98% cotton/2% elastane
- Latex-free

Item no. 45105 White, 5 cm x 4.5 m stretched

All fixation bandages freely combinable

Self-fixing



✓ MEDICAL PRODUCT CE

PINOFIT® cohesive bandage

- For support and fixation bandages
- Does not stick to skin or hair
- Also usable in combination with classic tape
- Elastic, self-adherent
- 98% cotton/2% elastane
- Latex-free

Item no. 45106 White, 5 cm x 4.5 m stretched



✓ MEDICAL PRODUCT CE

LEUKOTAPE® classic

- For classic tape bandages
- Inelastic adhesive bandage with zinc oxide rubber adhesive
- High adhesiveness
- Cotton fabric

Item no. 33580 12 rolls (3.75 cm x 10 m)

MASSAGE

PINO MASSAGE
THE THERAPIST'S CHOICE
SINCE 1922



SPORT

Ideal for soothing sports massages before or after training.

MADE
IN
GERMANY



DERMATOLOGICALLY TESTED



LABORATORY-TESTED

Massage lotion Sport Arnica

- Arnica flower extract to enable massaging as potential regeneration support after sports injuries

Item no. 35650 1000 ml bottle

Massage lotion Sport Magnesium

- With magnesium chloride for massage after exercise
- Specially designed for professional massage applications

Item no. 35630 1000 ml bottle


Fascia wax Sport

- Perfect for selective massages and intensive fascia treatments
- Provides especially intense nourishment of the skin with the addition of beneficial arnica extract, apricot kernel oil, almond oil and vitamin E

Item no. 32695 150 ml can

WHICH MASSAGE PRODUCT IS RIGHT FOR YOU?

Using our viscosity scale, you can easily determine the perfect massage product for your needs!



*Property increases in the direction of the arrow

POMEGRANATE ARGAN OIL

Fruity and fresh with optimal grip.

The fruity-fresh scent of Mediterranean pomegranates invigorates body and mind during the massage. Thanks to its hydrating ingredients, beneficial argan oil nourishes and cares for all skin types, leaving skin feeling supple and soft.



SENSITIVE

Especially suitable for sensitive skin and also for allergy sufferers.

The Sensitive range of massage products was especially developed for the needs of particularly sensitive skin. They are free from perfumes, preservatives and dyes and provide the skin with intensive moisture.



TOP
SELLER



Our highlights
for your massage

- ⊕ Natural argan oil with skin nourishing properties
- ⊕ Nourishes the skin during massage
- ⊕ Special massage pumps for dispensing on p. 93

Supplied without pump.



Massage lotion

- Extra light texture specially developed for all types of applications
- Perfect balance of grip and slipperiness

Item no. 35370 1000 ml bottle

Item no. 35375 10 L canister

Massage oil

- Extreme viscosity for maximum glide
- Specially developed for medium to large area massage applications

Item no. 35380 1000 ml bottle

Item no. 35385 10 L canister

Massage milk

- Fluid, slippery texture
- Specially developed for medium to large area massage applications

Item no. 35315 1000 ml bottle

Massage creme

- Low viscosity for high grip
- Ideal for small area applications

Item no. 35395 1000 ml can



Our highlights
for your massage

- ⊕ Free from preservatives, dyes and perfume
- ⊕ Suitable for allergy sufferers
- ⊕ Special massage pumps for dispensing on p. 93

Supplied without pump.



Massage oil

- Specially developed for large and elongated applications

Item no. 35530 1000 ml bottle

Item no. 35567 10 L canister

Massage ointment

- Replenishing wool wax is particularly suitable for sensitive skin
- Perfect for spot massages

Item no. 35414 500 ml can

Massage lotion

- Vitamin-rich avocado oil and beneficial shea butter
- Suitable for all types of applications

Item no. 35535 1000 ml bottle

Item no. 35538 10 L canister

Massage milk

- Skin-conditioning properties thanks to beneficial oat oil and antioxidant vitamin E
- Optimal for medium to large area applications

Item no. 35305 1000 ml bottle

Massage creme

- With natural shea butter extract
- High grip for small-area applications

Item no. 35540 1000 ml can



AVOCADO

Intensive skin care with vitamin-rich avocado oil.

The rich avocado oil intensively nourishes stressed skin. It provides a lot of moisture for dry skin in particular, without leaving a greasy film.

MADE
IN
GERMANY

ALOE VERA

With the care formula for velvety soft skin.

Aloe Vera is known for its moisturising and antioxidant properties. The velvety massage products with their creamy scent and aloe vera extract nourish the skin of therapists and patients alike.

MADE
IN
GERMANY



Our highlights
for your massage

- ⊕ Natural avocado oil with skin nourishing properties
- ⊕ Provides dry skin with plenty of moisture
- ⊕ Special massage pumps for dispensing on p. 93

Supplied without pump.



Our highlights
for your massage

- ⊕ Natural aloe vera extract with moisturising properties
- ⊕ Provides dry skin with plenty of moisture
- ⊕ Special massage pumps for dispensing on p. 93

Supplied without pump.

Massage creme

- Low viscosity for high grip
- Ideal for small area applications

Item no. 35316 1000 ml can

Massage milk

- High viscosity for good glide
- Optimal for medium to large area applications

Item no. 35365 1000 ml bottle

Massage lotion

- Perfect balance of grip and slipperiness
- Ideal for all types of application from spot to large area

Item no. 35156 1000 ml bottle

Item no. 35129 10 L canister

Massage oil

- Also contains pure peanut oil for maximum spreadability
- Perfect for large-area massages

Item no. 35160 1000 ml bottle

Item no. 35138 10 L canister

Massage creme

- Low viscosity for high grip
- Ideal for small area applications

Item no. 35781 1000 ml can

Massage milk

- High viscosity for good glide
- Especially light texture for medium to large area massages

Item no. 35345 1000 ml bottle

Massage lotion

- Especially beneficial for skin care thanks to aloe vera, panthenol (provitamin B5) and allantoin
- Optimal for all types of applications

Item no. 35192 1000 ml bottle

Item no. 35147 10 L canister

Massage oil

- Intensively nourishes the skin with aloe vera extract and sunflower oil
- Perfect for large-area massage applications

Item no. 35182 1000 ml bottle

Item no. 35137 10 L canister



LEMON

With the tangy fragrance of fresh lemons.

This fragrance invigorates the senses! An application with massage products from the Lemon series can have an activating and refreshing effect and not only features a fresh and tangy scent but also provides optimal care for the skin.

MADE
IN
GERMANY



Our highlights for your massage

- ⊕ Fresh tangy fragrance
- ⊕ Nourishes the skin during massage
- ⊕ Special massage pumps for dispensing on p. 93

Supplied without pump.

Massage creme

- Low viscosity for high grip
- Perfect for small area applications

Item no. 35780 1000 ml can

Massage milk

- With nourishing sunflower oil and antioxidant vitamin E
- Optimal for medium to large area massages

Item no. 35335 1000 ml bottle

Massage lotion

- With vitamin-rich avocado oil
- Optimal for all types of applications

Item no. 35285 1000 ml bottle

Item no. 35295 10 L canister

Massage oil

- Extreme viscosity for maximum glide
- Perfect for large-area massage applications

Item no. 35460 1000 ml bottle

Item no. 35405 10 L canister

DERMATOLOGICALLY TESTED

LABORATORY-TESTED

HEAT

With deep formula for comfortable, warm skin feeling.

Our heat massage products generate an intense, pleasant deep heat during the massage and are especially recommendable after sport activities.

MADE
IN
GERMANY



Massage ointment

- With natural chilli extract, which can support blood circulation
- Ideal for spot massages after sports

Item no. 35503 150 ml can

Item no. 35496 500 ml can

Please follow the instructions on the label.

DERMATOLOGICALLY TESTED

LABORATORY-TESTED

Massage lotion

- With stimulating rosemary extract and warming vanillyl butyl ether
- Suitable for all types of applications

Item no. 35245 1000 ml bottle



Our highlights for your massage

- ⊕ Intensive deep warmth
- ⊕ Ideal for massages after sports activities
- ⊕ Special massage pumps for dispensing on p. 93

PURE

Ideal for sensitive skin.

Massage lotion

- Especially suited for the treatment of sensitive skin
- Without perfume, without dyes, free from parabens and paraffin
- Moisturising ingredients help to keep the skin supple

Item no. 35162 1000 ml bottle



DERMATOLOGICALLY TESTED

LABORATORY-TESTED

SPORT AND MASSAGE CREAMS

Therapeutic sports massages are designed to improve performance, regeneration and muscle recovery in athletes. Depending on the desired massage effect, warming or cooling ointments can help to relax or stimulate targeted muscle groups and joints. Our range of massage ointments has been specially developed for different needs in therapy, training and competition.

PINOFIT® heat balm

With chilli. Increases blood circulation in skin and warms.

- With noticeable deep formula for comfortable, warm skin feeling
- The added beneficial calendula extract, soothing shea butter, nourishing sunflower oil and vitamin E contribute to the care of the skin

Item no. 32300 90 ml tube

PINOFIT® leg balm

For tired legs. Reactivating & relaxing.

- Perfect for the cooling massage after sport and to pamper tired and heavy legs
- Chestnut extract, well-known for its vein-strengthening properties, combined with red vine extract and mint oil is ideal for a pleasant pamper session

Item no. 32690 90 ml tube

PINOFIT® arnica cream

With arnica extract. Nourishing & invigorating.

- Contains high quality arnica extract and nourishes the skin, making it supple
- Carefully obtained Arnica montana extract stands out for its anti-inflammatory, revitalising and soothing properties
- Also recommended for a relaxing massage

Item no. 32400 90 ml tube

PINOFIT® cooling gel

With freshness effect. Revitalising & cooling.

- Cools pleasantly with mint oil and menthol
- The freshness formula nourishes and revitalises heavy legs and helps to relax and recover after exercises

Item no. 32200 90 ml tube

All 90 ml varieties freely combinable



Supplied without creams.

PINOFIT® display

- High-quality display for presenting PINOFIT® sport and massage creams
- Space-saving, for individual display configurations

Item no. 32540 Cardboard box, set up approx. 22.5 x 19.5 x 7.5 cm

MADE
IN
GERMANY

BODY TONIC

The ideal way to start or finish any massage. Refreshes the skin and removes residual massage oil.



All 1000 ml varieties freely combinable

Pink Grapefruit

With 70% vol. ethanol.

With the aromatic scent of pink grapefruit.

Item no. 35550 1000 ml bottle



Lemon Thyme

With 96% vol. ethanol.

With refreshing lemon & thyme fragrance.

Item no. 34570 1000 ml bottle



Orange Basil

With 96% vol. ethanol.

With Mediterranean orange & basil fragrance.

Item no. 34668 1000 ml bottle

Also available in 10-litre

Item no. 34677 10 L canister

MASSAGE BUTTER

An experience for the body and senses, with rich cocoa and shea butter.



All varieties freely combinable



COMPOSITION OF RICH COCOA & SHEA BUTTER AS WELL AS JOJOBA OIL AND VITAMIN E

100% NATURAL

Orange Spirit – balance –

The enigmatic scent of the essential oils of West Indian sandalwood and aromatic oranges dissipates a cosy warmth.

Item no. 34210 100 ml container

Cherry Blossom – calm –

The gentle scent of precious cherry blossom extract contributes to the balance between mind and body.

Item no. 34205 100 ml container

Orange Lemongrass – active –

The aromatic fragrance of orange oil and refreshing lemongrass oil gives the massage energising properties.

Item no. 34215 100 ml container

MADE
IN
GERMANY

NATURAL AROMA MASSAGE OILS

For the highest standards in a massage.




All 500 ml varieties freely combinable



All 3 L varieties freely combinable

MADE
IN
GERMANY

 HIGH-QUALITY CARE COMPOSITION OF SMOOTHING SOYBEAN OIL,
SOOTHING COCONUT OIL & PROTECTIVE JOJOBA OIL

 100% NATURAL

Lavender Mallow – calm –

Lavender and mallow are considered relaxing and stress-reducing and can support the balance of body and mind.

Item no. 38146 500 ml glass bottle with pump | *Item no. 38807* 3 L bag-in-box

Green Tea Ginger – anti cellulite –

Green tea and ginger are known for their stimulating and antioxidant properties and give the massage oil its stimulating character.

Item no. 38153 500 ml glass bottle with pump | *Item no. 38853* 3 L bag-in-box

Orange Lemongrass – active –

Activating thanks to its combination of fragrant orange and invigorating lemongrass oil.

Item no. 38110 500 ml glass bottle with pump | *Item no. 38810* 3 L bag-in-box

Lime Bamboo – active –

Revitalising lime scent and skin -protecting bamboo extract invigorate body and mind.

Item no. 38150 500 ml glass bottle with pump | *Item no. 38850* 3 L bag-in-box

Cherry Blossom – balance –

Delicately scented extracts of cherry blossom make skin radiant and give the massage oil its balancing character.

Item no. 38108 500 ml glass bottle with pump | *Item no. 38890* 3 L bag-in-box

Orange Spirit – balance –

The fragrance of fresh oranges and the enigmatic scent of West Indian sandalwood can support the balance between body and mind during the massage.

Item no. 38112 500 ml glass bottle with pump | *Item no. 38855* 3 L bag-in-box

Papaya Aloe – balance –

Balancing thanks to the exotic fruity fragrance of papaya and moisturising aloe vera extract.

Item no. 38143 500 ml glass bottle with pump | *Item no. 38867* 3 L bag-in-box

TOP
SELLER

Lemon Hibiscus – active –

The refreshing scent of lemon oil and precious extracts of hibiscus flower invigorate the senses.

Item no. 38103 500 ml glass bottle with pump | *Item no. 38865* 3 L bag-in-box

Swiss Pine – calm –

The aromatic-woody scent of pine oil has a calming effect and supports deep relaxation.

Item no. 38133 500 ml glass bottle with pump | *Item no. 38875* 3 L bag-in-box

Coconut Mint – balance –

With especially high nourishing coconut oil content and stimulating mint oil for balancing properties.

Item no. 38180 500 ml glass bottle with pump | *Item no. 38885* 3 L bag-in-box

Wild Rose – calm –

Exquisite rose petal extract and luxurious wild rose oil give this massage oil calming properties.

Item no. 38160 500 ml glass bottle with pump | *Item no. 38860* 3 L bag-in-box

Lemon Mint – active –

Refreshing lemon and stimulating mint oil give the massage oil its invigorating character.

Item no. 38100 500 ml glass bottle with pump | *Item no. 38893* 3 L bag-in-box

Mother-of-Pearl – skin protect –

With nourishing vital substances, remineralising mother-of-pearl extract and jojoba, coconut and soybean oils for supple skin.

Item no. 38675 500 ml glass bottle with pump

Honey Ginger – balance –

Gently stimulating ginger and skin-calming honey extract give this massage oil balancing properties.

Item no. 38170 500 ml glass bottle with pump | *Item no. 38870* 3 L bag-in-box



NATURAL AROMA MASSAGE CANDLES

For soothing warmth and care during massage. With nourishing coconut oil, sunflower oil and beeswax.

MADE
IN
GERMANY



All 60 g varieties freely combinable

 BEESWAX AND COCONUT OIL ARE INTENSIVE SKIN-CARE COMPONENTS

 100% NATURAL

Lavender Mallow – calm –

In addition to the pleasantly nourishing beeswax, the herbal-fresh scent of lavender and mallow contributes to holistic relaxation.

Item no. 35024 60 g can

Green Tea Ginger – anti cellulite –

Green tea and ginger are known for their stimulating and antioxidant properties and give the massage candle its stimulating character.

Item no. 35053 60 g can

Orange Lemongrass – active –

Caresses with the stimulating fragrance of aromatic orange and lemongrass oil.

Item no. 35042 60 g can

Lime Bamboo – active –

Revitalising lime oil and skin-protecting bamboo extract invigorate body and mind.

Item no. 35052 60 g can

Cherry Blossom – balance –

Extracts of cherry blossom and skin-nourishing waxes and coconut oil provide a balancing character.

Item no. 35004 60 g can

Lemon Hibiscus – active –

The refreshing scent of lemon oil and precious extracts of hibiscus flowers invigorate the senses.

Item no. 35082 60 g can

Swiss Pine – calm –

The aromatic-woody scent of pine oil has a calming effect and supports deep relaxation.

Item no. 35092 60 g can

Coconut Mint – balance –

Beneficial coconut oil and stimulating mint oil nourish the skin and give the massage candle balancing properties.

Item no. 35032 60 g can

Wild Rose – calm –

Caresses with the gently soothing fragrance of precious rose petal extract and wild rose oil.

Item no. 35012 60 g can

Honey Ginger – balance –

Skin-soothing honey and gently stimulating ginger extract help balance and calm stressed skin.

Item no. 35062 60 g can

TOP
SELLER

TOP
SELLER



Massage candle set with 12 candles, 60 g each, and display

- Two candles each of types Honey Ginger, Coconut Mint, Orange Lemongrass, Lavender Mallow and Lime Bamboo
- Includes display

Item no. 35100



Thanks to the optimal melting temperature of 50–60 degrees, the liquid wax is ready to be used for a warming massage directly: simply drizzle over the hand onto the body and rub. The pleasant warmth, the gratifying wax and enchantingly fragrant essential oils make the candle wax massage a feel-good experience of a special kind!

Carsten Boss, Senior Therapist



TRADITIONAL MASSAGE OILS

High-quality oils with natural properties.

MADE
IN
GERMANY



DERMATOLOGICALLY
TESTED "VERY GOOD"

100% VEGAN

100% NATURAL



TOP
SELLER

Arnica – nurturing –

- Gently cares for particularly sensitive skin with valuable extract from arnica flowers
- Arnica is said to have analgesic and anti-inflammatory effects and is useful in the treatment of sports injuries
- The Arnica montana extract stands out for its revitalising and harmonising properties
- Suitable for any skin type, even for sensitive and dry skin

Item no. 38780 500 ml glass bottle with pump

Sesame – protecting –

- Made from 100% pure sesame oil
- Traditionally used for ayurveda massages
- Sesame oil is said to encourage blood circulation and have warming and cleansing properties
- Particularly suitable for mixed skin with a tendency to dryness as well as dry and mature skin

Item no. 38720 500 ml glass bottle with pump

Item no. 38725 3 L bag-in-box with dispenser

Shea Aloe – sensitive –

- Mild composition of shea butter and aloe vera
- With the nourishing blend of smoothing soy, soothing coconut and protective jojoba oil
- Nourishes very sensitive skin
- Supports the suppleness of the skin and intensively moisturises it

Item no. 35605 500 ml glass bottle with pump

Item no. 35606 3 L bag-in-box with dispenser

Almond – moisturising –

- Made from 100% pure almond oil
- Almond oil is well absorbed into the skin
- It is soothing and calms the skin
- It is particularly suitable for chapped, sensitive and dry skin of any age
- Classic massage oil for any skin type

Item no. 38710 500 ml glass bottle with pump

Jojoba – repairing –

- Made from 100% pure jojoba oil
- High quality massage oil from the seeds of the jojoba bush
- Jojoba oil is very well absorbed into the skin, leaving it smooth and protecting it from dehydration
- Also for sensitive and brittle skin

Item no. 38730 500 ml glass bottle with pump

St. John's wort – relaxing –

- Made from pure St John's Wort oil on an olive oil basis
- Traditionally used for Breussmassages
- St John's Wort oil relaxes and warms, making it perfect for back massages

Item no. 38700 500 ml glass bottle with pump



OIL BLENDS

Perfect for large-area massages.



Herbal stamp massages have been practised for over 2000 years! The unique combination of herbs, heat, nourishing oils and an intensive massage produces the greatest possible relaxation of body and mind. Since large area massages use a lot of oil, I recommend our speciality oil blends in large units – for a perfect price-to-massage ratio!”

Carsten Boss, Senior Therapist

MADE
IN
GERMANY



Odourless even when heated

Special oil blend for herbal stamps and hot stones

- Made from almond, coconut, nourishing calendula and sunflower oils
- Well suited for an intensive hot stones & herbal stamps massage
- Odourless & heat-resistant

Item no. 38750 3 L bag-in-box with dispenser

100% NATURAL



Special massage oil almond

- From pure almond oil on soya oil basis
- With high-quality almond oil and vitamin E
- Also for sensitive skin
- Perfect for large-area massages

Item no. 38770 5 L canister

100% NATURAL

Special pump for canister

Item no. 39967

HERBAL STAMP

Warm massages for muscle relaxation and relaxation.

Professionally bound



- Bound by hand
- Various purely natural fillings
- Resealable aroma freshness pack
- Material from 100% unbleached cotton

Herbal stamp herbs

Item no. 32475 1 pack of 2 stamps, each 150 g

Herbal stamp lemongrass

Item no. 32453 1 pack of 2 stamps, each 150 g

All 150 g varieties freely combinable



Herbal stamp coconut & rice

Item no. 32480 1 pack of 2 stamps, each 200 g

Herbal stamp rose petals & rice

Item no. 32483 1 pack of 2 stamps, each 200 g



All 200 g varieties freely combinable

DISPENSER BOTTLES FOR MASSAGE PRODUCTS

For the quick use of small amounts of massage oil.



Glass dispenser bottle with lid and drip guard

Item no. 30410 170 ml



Plastic dispenser bottle with flip top cap

Item no. 39761 200 ml



Plastic dispenser bottle with pump

Item no. 39766 200 ml



SKIN CARE



SHOWER ME! SHOWER FOAMS

Creamy care foam with airy volume and the most beautiful scents of nature.

MADE
IN
GERMANY



Our highlights
for your skin

- ⊕ Rich & concentrated
- ⊕ Very creamy nourishing foam with full and airy volume



All 200 ml varieties freely combinable



DERMATOLOGICALLY TESTED

Pink Grapefruit

- With the Mediterranean scent of pink grapefruit
- The creamy foam with vitamin-rich sunflower oil leaves your skin feeling wonderfully conditioned

Item no. 37478 200 ml can

Cherry Blossom

- With the indulgent scent of cherry blossom
- Vitamin-rich sunflower oil luxuriously nourishes the skin and protects it from drying out

Item no. 37440 200 ml can

Lemon Tonic

- With the refreshing fragrance of tonic water with a splash of lemon
- Vitamin-rich sunflower oil luxuriously nourishes the skin and protects it from drying out

Item no. 37470 200 ml can

Lemon Hibiscus

- With the refreshing scent of lemons and hibiscus blossoms
- The creamy foam with vitamin-rich sunflower oil leaves your skin feeling wonderfully conditioned

Item no. 37432 200 ml can

Aqua & Minerals

- With the reviving fragrance of a soft ocean breeze
- The creamy foam with vitamin-rich sunflower oil leaves your skin feeling wonderfully conditioned

Item no. 37438 200 ml can

Passion Fruit

- With the exotic scent of delicious passion fruit
- Vitamin-rich sunflower oil pampers and gently nourishes the skin

Item no. 37441 200 ml can

Orange Olive Oil

- With the fragrance of sun-kissed oranges
- Exquisite olive oil nourishes the skin and makes it feel like velvet

Item no. 37480 200 ml can

Sea Algae and Water Mint

- With the touch of a fresh sea breeze
- The creamy foam with vitamin-rich sunflower oil leaves your skin feeling wonderfully conditioned

Item no. 37430 200 ml can

Mango Macadamia

- With the intense fragrance of ripe mango
- Exquisite macadamia oil luxuriously nourishes the skin and protects it from drying out

Item no. 37475 200 ml can

Pomegranate Argan Oil

- With the reviving fragrance of Mediterranean pomegranates
- Argan oil luxuriously nourishes the skin and protects it from drying out

Item no. 37485 200 ml can

TOP
SELLER



Shower Me! Shower foam advertising display

- Mobile advertising stand for placement in physiotherapy practice or spa
- Suitable for the presentation of a SHOWER ME! Shower foam of your choice (200 ml can)
- Can be configured to suit your needs

Item no. 37170 Cardboard box, set up: approx. 23 x 25 x 10 cm

Supplied without contents.

PROTECT ME! HAND CREAM CONCENTRATE

MADE
IN
GERMANY



Our highlights for your hands

- ⊕ With a particularly high shea butter content of more than 15%
- ⊕ Thanks to a special care formula, the hand cream is absorbed instantly, leaving your hands wonderfully smooth and soft



DERMATOLOGICALLY TESTED "VERY GOOD"

Pomegranate Argan Oil

- With a particularly high shea butter content of more than 15%
- With the Mediterranean fragrance of fresh pomegranates
- The special care formula provides silky soft hands

Item no. 36014 50 ml tube

Shea Butter

- With a particularly high shea butter content of more than 15%
- The innovative care formula with olive oil, panthenol and hyaluronic acid offers highly concentrated protection for very dry and stressed hands

Item no. 36016 50 ml tube



All varieties of hand cream and foam are freely combinable

HAND CREAM CARE FOAM



Our highlights for your hands

- ⊕ The unique Protect Me! Formula pampers the skin and nourishes hands silky smooth
- ⊕ Very creamy nourishing foam with full and airy volume
- ⊕ Rich and concentrated



DERMATOLOGICALLY TESTED



100% VEGAN

Aronia Berry

- Protects hands with macadamia nut oil, shea butter and urea
- Nourishes the skin with natural antioxidants of aronia berry extract

Item no. 33325 50 ml can

Pomegranate Argan Oil

- Protects hands with argan oil and shea butter
- With the Mediterranean fragrance of fresh pomegranates

Item no. 33335 50 ml can

Cherry Blossom

- Protects hands with almond oil, shea butter and panthenol
- With the scent of delicate cherry blossoms

Item no. 33310 50 ml can



PROTECT ME! BODY LOTIONS

For smooth and soft skin and long-lasting moisturisation.

MADE
IN
GERMANY



All varieties freely combinable



DERMATOLOGICALLY TESTED



100% VEGAN



HIGH QUALITY INGREDIENTS

Cherry Blossom

- Velvety care experience for the skin with the delicate scent of cherry blossoms
- The special PROTECT ME! Formula with shea butter, cherry blossom extract and hyaluronic acid cares for and intensively moisturises dry and stressed skin

Item no. 33338 200 ml tube

Pomegranate Argan Oil

- Special skin care with the Mediterranean scent of ripe pomegranates
- The special PROTECT ME! Formula of beneficial argan oil, rich shea butter and moisturising hyaluronic acid leaves skin feeling optimally cared for and velvety soft

Item no. 33330 200 ml tube

BALZAMON BEESWAX BALM

Rich skincare series based on beeswax with high quality ingredients. Handmade in the south of France.

All varieties freely combinable



100% NATURAL



Natural

Item no. 30355 30 ml glass jar

Pink

Item no. 30361 30 ml glass jar

Lavandula

Item no. 30356 30 ml glass jar

Hypericum

Item no. 30358 30 ml glass jar

Propolis

Item no. 30359 30 ml glass jar

- 100% natural with all ingredients from organic farming
- Without addition of fragrances, emulsifiers or preservatives
- During production and bottling by hand in the south of France, only direct sunlight is used for heating
- Suitable for skin care, rubs and massages
- High grip for spot application
- Massage with Balzamon can create a slightly warming sensation



NATURAL BODY OILS

Exclusive fragrance compositions with rich jojoba, coconut and soybean oil.

MADE
IN
GERMANY



Our highlights for your skin

- ⊕ With valuable extracts and essential oils
- ⊕ Suitable for all skin types, especially for dry and sensitive skin
- ⊕ Rich and concentrated

TOP
SELLER

Orange Lemongrass – active –

The blend of aromatic orange oil and refreshing lemongrass oil gives this unique massage oil energising properties.

Item no. 38615 100 ml bottle with pump



All 100 ml varieties freely combinable

⊕ HIGH-QUALITY CARE COMPOSITION OF SMOOTHING SOYBEAN OIL, SOOTHING COCONUT OIL & PROTECTIVE JOJOBA OIL

100% NATURAL

Lavender Mallow – calm –

Valuable essential oils unfold their unique aroma and give the body care its balancing character.

Item no. 38618 100 ml bottle

Lemon Hibiscus – active –

The blend of refreshing lemon oil and the fragrant hibiscus flower gives this unique body oil energising properties.

Item no. 38603 100 ml bottle

Cherry Blossom – balance –

Precious cherry blossom extracts give this pampering body oil balancing properties.

Item no. 38608 100 ml bottle

Coconut Mint – balance –

The blend of a high content in nourishing coconut oil and stimulating mint oil gives this unique body oil balancing properties.

Item no. 38610 100 ml bottle

Papaya Aloe – balance –

The blend of exotic fruity papaya fragrance and moisturising aloe vera extract gives this unique body oil harmonising properties.

Item no. 38643 100 ml bottle

Wild Rose – calm –

The blend of luxurious rose petal extract and exquisite wild rose oil gives this unique body oil calming properties.

Item no. 38600 100 ml bottle

Green Tea Ginger – anti cellulite –

Valuable ingredients nourish the skin and support it with an anti-cellulite effect.

Item no. 38653 100 ml bottle

Orange Spirit – balance –

The composition of essential oils of West Indian sandalwood and aromatic orange gives this unique body oil a balancing character.

Item no. 38612 100 ml bottle

Swiss pine – calm –

The aromatic woody scent of Swiss pine gives this soothing body oil calming properties.

Item no. 38633 100 ml bottle

Lime Bamboo – active –

The blend of revitalising lime fragrance and skin-protecting bamboo extract gives this unique body oil energising properties.

Item no. 38603 100 ml bottle

Honey Ginger – balance –

The blend of gently stimulating ginger and skin-calming honey extract gives this unique body oil balancing properties.

Item no. 38605 100 ml bottle

Mother-of-Pearl – skin protect –

The blend of nourishing vital substances and remineralising mother-of-pearl extract supports the smoothness of the skin and intensively rehydrates it.

Item no. 38670 100 ml bottle

TOP
SELLER

TOP
SELLER

BATH OILS

Pure relaxation with the forces of nature.

MADE
IN
GERMANY



All 1000 ml varieties freely combinable

All 10 L varieties freely combinable

 DERMATOLOGICALLY TESTED

 100% VEGAN

Pomegranate Argan Oil foam bath

- The scented foam bath with the delicious fragrance of fruity pomegranate relaxes and regenerates body and mind.
- The long-lasting body foam and exquisite argan oil provide intense skin care

Item no. 36056 1000 ml bottle

Bath oil Pine Needle

- High quality pine needle oil pleasantly supports muscle relaxation

Item no. 30599 1000 ml bottle

Breathe Freely bath

- Beneficial eucalyptus oil, menthol and pine needle essences promote easier breathing

Item no. 36082 1000 ml bottle

Bath oil Lavender

- Exquisite lavender oil pleasantly supports mental and physical relaxation

Item no. 31350 1000 ml bottle

Nurturing oil bath Almond Blossom

- Supports dry and sensitive skin with a high percentage of natural almond oil, leaving the skin feeling silky smooth

Item no. 31435 1000 ml bottle

Bath oil Balm

- Promotes soothing relaxation
- Exquisite Indian Melissa with its pleasant fragrance nourishes the skin and can offset negative tension

Item no. 31489 1000 ml bottle

Item no. 31498 10 L canister

Hay Flower Herbs

- Exquisite essences of hay flowers and meadow herbs pleasantly support muscle relaxation

Item no. 31185 1000 ml bottle

Item no. 31676 10 L canister



PRO SAUNA CONCENTRATES

Exclusive infusions with natural essential oils.

MADE
IN
GERMANY



Our highlights for your sauna

- ⊕ Suitable for sauna and steam bath
- ⊕ Evaporates without residues
- ⊕ For up to 100 infusions



All varieties freely combinable

Orange Lemongrass

Activates body and soul with aromatic orange and reviving lemongrass oil.

Item no. 37510 1000 ml bottle

Lemon Mint

Revitalises the body with refreshing lemon oil and gently stimulating mint oil.

Item no. 37545 1000 ml bottle



Rosemary Eucalyptus

Encourages the regeneration of body and soul with refreshing eucalyptus and warming rosemary oil.

Item no. 37580 1000 ml bottle

Swiss Pine

Supports deep relaxation of body and soul with relaxing, aromatic, woody Swiss pine.

Item no. 37592 1000 ml bottle

Eucalyptus Menthol

Revives body and soul with refreshing eucalyptus oils and menthol.

Item no. 37338 1000 ml bottle



THERMOTHERAPY

MUD PACKS

Reliable top quality for professional applications.



LABORATORY-TESTED



100% NATURAL



MEDICAL PRODUCT



PINO natural mud pack

- Ready-to-use upland moor mud pack with balanced moisture content and a high content of substrate
- Suitable for a localised thermotherapy
- Processing is 100% additive free; the valuable natural humic acid remains unchanged
- The high content of dry matter makes it an excellent hygroscopic agent, facilitates a very high tolerance of fluctuations in humidity and guarantees an ideal consistency for therapy
- Permeable special gauze for intensive skin contact
- With a useful perforation to separate moor mud, gauze and foil for an environment friendly disposal
- Clean and simple application
- Recommended in combination with PINO heat carriers
- Disposable pack

Pack of

Moor mud pack size 1

Item no. 30085 40 packs, each 38 x 28 cm

Pack of

Moor mud pack size 2

Item no. 30080 20 packs, each 58 x 38 cm

Pack of

Moor mud pack sz. 1, PREMIUM, extra thick

Item no. 30075 30 packs, each 38 x 28 cm

Our highlights for your application

- ⊕ With valuable natural humic acids
- ⊕ Ready for use with practical perforation for easy separation of moor mud, gauze & foil
- ⊕ Permeable special gauze for intensive skin contact

HEAT CARRIERS

Reliable top quality for professional applications.



Our highlights for your application

- ⊕ Can be perfectly modelled for any application
- ⊕ Ideally suited to heating PINO moor mud, hay, Swiss pine and bee wax packs
- ⊕ Impresses because of its high thermal storage capacity and extremely consistent heat emission
- ⊕ Simple to use and clean

Water baths see p. 125

✓ MEDICAL PRODUCT CE



Heat carrier with moor mud filling


- Impresses because of its high thermal storage capacity and extremely consistent heat emission
- Due to naturally moist moor filling optimum mouldability on all treatable body parts
- Simple to use and clean
- To be heated in a water bath
- Not suitable for heating cabinets or microwave

Size 1

Item no. 30470 For mud pack size 1 (27.5 cm x 37 cm)

Size 2

Item no. 30475 For mud pack size 2 (57 cm x 37 cm)



MADE
IN
GERMANY

Heat carrier with paraffin filling

- The melting point of paraffin ensures a long lasting heat distribution in the area to be treated during the whole lengths of the treatment
- Incorporated gauze ensures the equal distribution of the filling
- Highly robust due to the high quality and extremely flexible foil material
- Simple to use and clean
- To be heated in a water bath or heating cabinet
- Not suitable for microwave

Size 1

Item no. 30374 For mud pack size 1 (27.5 cm x 37 cm)

Size 2

Item no. 30376 For mud pack size 2 (57 cm x 37 cm)



MADE
IN
GERMANY

Heat carrier with special fluid

- Extraordinarily light and thin heat carrier
- Flexible and malleable even when cold
- Incorporated gauze ensures the equal distribution of the filling
- Simple to use and clean
- To be heated in a water bath
- Not suitable for heating cabinets or microwave

Size 1

Item no. 30375 For mud pack size 1 (27.5 cm x 37 cm)


Size 2

Item no. 30017 For mud pack size 2 (57 cm x 37 cm)


PINOTHERM® SPECIAL HEAT PACKS

Energy-saving, as no fango blocks need to be melted.

MADE IN GERMANY




Special chamber system



MEDICAL PRODUCT CE

PINOTHERM® THERMAL ENERGY PROGRESS



Long lasting heat emission

- Filled with paraffin-coated PCM silicate (phase change material)
- Particularly long-lasting and even heat emission
- Fully compensable by all state health insurance providers in accordance with HMR s. 24 (2) No. 5

Practical in use

- Smooth surface without any irritating seams
- Time-saving and easy to heat
- Ideally suited for the warming cabinet
- Perfect application and malleability properties
- Easy to maintain and clean
- Easy to disinfect

Professional in quality

- High heat storage capacity
- Stable double stitching of the chamber system
- Material made of skin-friendly and durable nylon fabric

Smooth lying surface for more comfort

- Optimised, flat surface means greater comfort for your patients
- Facilitates residue-free disinfection during wiping and thus protects the material
- System of special chambers prevents slipping of the silicate and thus ensures uniform heat distribution

- | | | |
|--|---|--|
| Small
Item no. 30250 1 pack approx. 30 x 50 cm | Medium
Item no. 30270 1 pack approx. 40 x 60 cm | Large
Item no. 30260 1 pack approx. 50 x 70 cm |
|--|---|--|

FANGO

Heat therapy with fango in proven top quality.

Fango is a volcanic mineral mud that is very good at storing and transferring heat. Heat therapy with fango causes relaxation of muscles and promotes blood circulation. PINOTHERM fangoparaffin is also mixed with kerosenes with different melting points. This makes it easy to model and particularly uncomplicated to use. Try it out!

Carsten Boss, Senior Therapist



MEDICAL PRODUCT CE

PINOTHERM® fangoparaffin

- With the characteristics of prehistoric fossilised sea mud (Posidonia Shale), it has fossil content such as minerals and petroleum
- Excellent heat storage and heat transfer capacity
- Fast melting & easy to model

The two varieties are freely combinable

Fangoparaffin type N for low viscosity
Item no. 30400 24 bars each 500 g = 12 kg

Fangoparaffin type K for higher viscosity
Item no. 30419 24 bars each 500 g = 12 kg

Anti-static

Disposable fango foil rolled

- Temperature-resistant up to 100°C

Item no. 30142 15 µ thick, roll 1000 m x 65 cm

Item no. 30160 25 µ thick, roll 750 m x 65 cm

Anti-static

Disposable fango foil on a tear-off pad

Item no. 30124 25 µ thick, 10 x 100 sheets, 55 x 75 cm

Item no. 30151 25 µ thick, 10 x 100 sheets, 55 x 85 cm

Fango cloth (not shown)

- 100% cotton
- Suitable for tumble dryer up to 60°C

Set of 2

Item no. 44879 2 pcs., 250 x 150 cm, approx. 1650 g, white

Fango sheet (not shown)

- 100% cotton
- Suitable for tumble dryer up to 60°C

Set of 2

Item no. 44860 2 pcs., 250 x 150 cm, approx. 1150 g, white

Disposable gauze (not shown)

Item no. 40460 100 m, width 65 cm, white

Paraffin cleaner (not shown)

Item no. 39645 5 L canister

MADE IN GERMANY

WELLNESS BODY WRAPS

Soothing warmth and aromatic scents.



100% NATURAL



100% NATURAL

Pack of 10

Swiss Pine wrap

- Filled with Swiss pine shavings for a soothing heat treatment
- The heat treatment with the Swiss pine wrap is part of the unique Swiss pine massage for intensive rejuvenation of body and mind
- The fragrant woody scent of the Swiss pine wrap stimulates while providing deep relaxation

Item no. 30091 10 wraps each 38 x 28 cm

Pack of 50

Hay wraps

- Ready to use hay wrap with fermented, natural hay
- For relaxation before or after the massage
- With a useful perforation to separate hay, gauze and foil for an environment friendly disposal.

Item no. 30081 50 wraps each 38 x 28 cm

100% NATURAL

Pack of 20

Beeswax wrap

- Ready to use beeswax wrap for honey massages
- Provides a lasting velvety soft, supple skin feeling
- Added beneficial antimicrobial propolis revitalises stressed skin

Item no. 30350 20 disposable wraps each 45 x 33 cm

MADE IN GERMANY



HEAT PACK WITH VOLCANIC GRANULES

Top quality Made in Germany.



Four chambers provide ideal support and distribution of the granules

MADE IN A NON-PROFIT INCLUSION COMPANY

For long-lasting deep heat



Volcanic Pack

- Professional heat pack, also suitable for home use
- With natural volcanic granules for long-lasting warmth
- Simple and fast to use
- Heat in the microwave or oven
- Multiple reusability
- Handmade from Germany



Perfectly fits the shoulder area, with additional high area for the neck.

VULKAN PAK

New

Volcanic Pack

Item no. 30367 35 x 22 cm, 1045 g

Volcanic Pack Neck

Item no. 30368 42 x 31 cm, 2017 g

MEDICAL PRODUCT CE

PINOFIT® Cold Pack

- Offers optimum plasticity on all treatable body parts
- For cold therapy
- Ideal gel consistency
- Extremely soft and smooth, multi-layered special polyurethane foil for optimum flexibility
- Easy to clean with soap and water or with the alcohol-free SEPTA-PIN® AQUAQUICK rapid disinfection p. 104

Small

Item no. 30960 Approx. 13 x 28 cm

Medium

Item no. 30964 20 x 19 cm

Large

Item no. 30973 40 x 21 cm



MADE IN GERMANY

PRACTICE HYGIENE



NOSE SLIT TOWELS

The hygiene standard – suitable for any headrest.

Pack of 100

TOP
SELLER

Extra large disposable nose slit towels


- Clean, hygienic solution for your patients
- Protects the therapy bed from cosmetic products thus prolonging the life of the upholstery

Paper


Item no. 40850 1 pack, 100 pieces, extra large, approx. 40 x 30 cm

Extra soft fleece

Item no. 40865 1 pack, 100 pieces, extra large, approx. 40 x 30 cm



Paper



Fleece

All varieties freely combinable

Pack

Terry cloth nose slit towels

- Very fluffy Terry cloth quality approx. 500 g/m²
- 100% cotton
- Clean and hygienic cover
- Washable up to 60°C
- Suitable for tumble dryers
- Tested for harmful substances according to OEKO-TEX® Standard 100



Lime

Item no. 40436 Extra large, 50 cm x 60 cm

Grey

Item no. 40438 Extra large, 50 cm x 60 cm

Saffron

Item no. 40437 Extra large, 50 cm x 60 cm

Light grey

Item no. 40432 Extra large, 50 cm x 60 cm



Wall mount fixture

- Suitable for nose slit towels Approx. 40 x 30 cm
- Practical and clean
- Powder-coated sheet steel
- Supplied without fixtures and products

Silver

Item no. 48195 Approx. 51.5 x 41 x 16 cm

SPECIAL PUMPS & DISPENSERS

Clean filling and easy dosing.



Pump for massage cream container

- Suitable for all 1000 ml cans of massage products
- Adjustable dosage volume
- Especially clean thanks to self-closing protective cap
- The dispenser can be closed by turning the pump head

Item no. 39980

Pump for massage lotion bottles

Suitable for all 1000 ml bottles of massage products

Item no. 39216

Pump for canister

- Suitable for all 10 L canisters of massage products, SEPTAPIN® wash lotions, bath oil and body tonic
- Suitable for 5 L canister of special massage oil Almond

Item no. 39967

FOR PORTABLE MASSAGE BENCHES AND CHAIRS

With elastic band for perfect fit.



- Disposable round covers from white fleece
- Clean and hygienic cover
- Protects upholstery from cosmetic products
- Air- and water-permeable
- The elastic round cover fits the shape of the headrest perfectly and stays in place

Item no. 44275 1 pack, 100 pieces

Pack of 10



Disposable sheets for therapy beds from white fleece

- Clean and hygienic cover
- Protects upholstery from cosmetic products
- Air- and water-permeable
- Secured fit with elastic

Item no. 44265 1 pack, 10 pieces, approx. 210 x 80 cm

Pack of 50



Fleece sheets, washable

- Made from high quality fleece
- Can be washed up to 200 times up to 60°C
- Tear-resistant, lint-free, non-iron
- Ideal for fango treatments
- For all PINO therapy beds and portable massage benches

Item no. 40465 1 pack, 50 pieces, approx. 200 x 80 cm

PAPER TOWELS

Cleanliness & hygiene for you & your patients.



Paper towels, folded, natural white

- Soft and absorbent
- 2-ply quality
- Skin-friendly
- V-fold

Natural white
25 x 23 cm, 2-ply, V-fold,
bulk pack of 3200

Item no. 39785



Folded hand towels, blue-white serrated, extra soft

- Soft and absorbent
- 2-ply quality
- Skin-friendly
- V-fold

White-blue serrated
24 x 21 cm, 2-ply, V-fold,
bulk pack of 3750

Item no. 39755



Paper towels, folded ECO Natural

- Soft and absorbent
- 2-ply quality, FSC®-recycled
- Undyed and skin-friendly
- V-fold

Natural
Approx. 25.3 x 21 cm, 2-ply, VZZ-fold
Bulk pack with 3800 pieces

Item no. 39768

TOWEL DISPENSER

For hygienic hand drying.

Towel rolls

- Soft and absorbent
- Skin-friendly
- Towel dispensation from inside

Bright white

Item no. 39795 Large pack with 12 rolls, 1-ply, 120 metres per roll



Roller towel dispenser

- Sturdy, grey-white plastic case
- With window on the side for checking filling level
- Lockable
- For mounting on wall

White-grey

Item no. 39790 Plastic, 184 x 336 x 167 mm



Dispenser for folded towels

- Compact, space-saving form
- Sturdy grey-white plastic housing with side viewing window to check fill level
- Lockable
- For mounting on wall
- Room for roughly 320 sheets

Standard for 320 sheets

Item no. 39700 Plastic white-grey, 290 x 290 x 160 mm

Large for 480 sheets

Item no. 39800 Plastic white-grey, 450 x 290 x 160 mm



You can find mouth guards and disposable gloves

PRO COTTON FITTED SHEETS

Extra soft, comfortable & durable.

Extra strong professional quality –



Professional fitted sheets from fluffy terry cloth

- Especially soft
- Highly hard-wearing and durable
- Excellent fit
- Surface made of pure cotton
- 80% cotton and 20% polyester
- Quality: approx. 220 g/m²
- Suitable for tumble dryers
- Hem stitched elastic
- Washable up to 60°C
- Tested for harmful substances according to OEKO-



65 cm wide without nose slit

Saffron <i>Item no. 40401</i> Approx. 195 x 65 cm	Light beige <i>Item no. 40479</i> Approx. 195 x 65 cm
White <i>Item no. 40402</i> Approx. 195 x 65 cm	Light grey <i>Item no. 40408</i> Approx. 195 x 65 cm
Lime <i>Item no. 40407</i> Approx. 195 x 65 cm	Dark grey <i>Item no. 40413</i> Approx. 195 x 65 cm



Optimized with double stitching

65 cm wide with nose slit


Saffron <i>Item no. 40501</i> Approx. 195 x 65 cm	Light beige <i>Item no. 40508</i> Approx. 195 x 65 cm
White <i>Item no. 40502</i> Approx. 195 x 65 cm	Light grey <i>Item no. 40548</i> Approx. 195 x 65 cm
Lime <i>Item no. 40547</i> Approx. 195 x 65 cm	Dark grey <i>Item no. 40514</i> Approx. 195 x 65 cm

80 cm wide without nose slit

Saffron <i>Item no. 40489</i> Approx. 195 x 80 cm	Light beige <i>Item no. 40497</i> Approx. 195 x 80 cm
White <i>Item no. 40492</i> Approx. 195 x 80 cm	Light grey <i>Item no. 40428</i> Approx. 195 x 80 cm
Lime <i>Item no. 40427</i> Approx. 195 x 80 cm	Dark grey <i>Item no. 40416</i> Approx. 195 x 80 cm

80 cm wide with nose slit

Saffron <i>Item no. 40518</i> Approx. 195 x 80 cm	Light beige <i>Item no. 40528</i> Approx. 195 x 80 cm
White <i>Item no. 40519</i> Approx. 195 x 80 cm	Light grey <i>Item no. 40538</i> Approx. 195 x 80 cm
Lime <i>Item no. 40537</i> Approx. 195 x 80 cm	Dark grey <i>Item no. 40516</i> Approx. 195 x 80 cm



SHORT FITTED SHEETS

Without headrest.

Professional fitted sheets SHORT

- Very soft
- Highly hard-wearing and durable
- Excellent fit
- Extra strong professional quality
- Surface made of pure cotton
- 80% cotton and 20% polyester
- Quality: approx. 220 g/m²
- Suitable for tumble dryers
- Hem stitched elastic
- Washable up to 60°C
- Tested for harmful substances according to OEKO-TEX® Standard



65 cm wide

Lime <i>Item no. 40550</i> Approx. 160 x 65 cm	Saffron <i>Item no. 40552</i> Approx. 160 x 65 cm
Dark grey <i>Item no. 40551</i> Approx. 160 x 65 cm	Light grey <i>Item no. 40553</i> Approx. 160 x 65 cm

80 cm wide

Lime <i>Item no. 40556</i> Approx. 160 x 80 cm	Saffron <i>Item no. 40558</i> Approx. 160 x 80 cm
Dark grey <i>Item no. 40557</i> Approx. 160 x 80 cm	Light grey <i>Item no. 40554</i> Approx. 160 x 80 cm



COTTON HAND TOWELS





Highest professional quality for your practice.



For perfect cleanliness during your treatments

Top-quality terry hand towels

- Very fluffy quality approx. 500 g/m²
- Carefully tested for harmful substances according to OEKO-TEX® Standard 100
- 100% cotton
- Highly absorbent, skin-friendly and durable
- Particularly durable and washable up to 60°C
- Suitable for tumble dryers
- Available as guest towels, hand towels, shower towels or sauna



Pack of

30 x 50 cm guest towels

Made of fleece terry.

Yellow

Item no. 37991 30 x 50 cm

White

Item no. 38001 30 x 50 cm

Grey

Item no. 38005 30 x 50 cm

Lime

Item no. 38003 30 x 50 cm



Royal blue

Item no. 38002 30 x 50 cm

Pack of 4

50 x 100 cm hand towels

Made of fleece terry.

Yellow

Item no. 37992 50 x 100 cm

White

Item no. 38006 50 x 100 cm

Grey

Item no. 38009 50 x 100 cm

Lime

Item no. 38008 50 x 100 cm



Royal blue

Item no. 38007 50 x 100 cm

Pack of 2

70 x 140 cm shower towels

Made of fleece terry.

Yellow

Item no. 37993 70 x 140 cm

White

Item no. 38011 70 x 140 cm

Grey

Item no. 38016 70 x 140 cm

Lime

Item no. 38014 70 x 140 cm



Royal blue

Item no. 38012 70 x 140 cm

Pack of 2

80 x 200 cm sauna towels

Made of fleece terry.

Yellow

Item no. 37994 80 x 200 cm

White

Item no. 38017 80 x 200 cm

Grey

Item no. 38021 80 x 200 cm

Lime

Item no. 38019 80 x 200 cm



Royal blue

Item no. 38018 80 x 200 cm

EQUIPMENT-ASSISTED THERAPY




PULLEY APPARATUS

For effective strength & coordination training.

MADE
IN
GERMANY



 You can find further accessories online.

✓ MEDICAL PRODUCT CE

PINOFIT® pulley apparatus

- Standard with two handles
- With various options for wall mounting
- Load: 4 x 0.5 kg, 11 x 2.0 kg, total weight 24 kg
- All pulleys are fitted with ball bearings
- The rope length can be adjusted via a serial perforation irrespective of the rope height
- Frame colour: silver-anthracite

Item no. 48120 W x D x H: 31 x 34 x 232 cm

PINOFIT® pulley apparatus Explosive with 2 handles

- Especially for fast weight training
- With various options for wall mounting
- Load: 10 x 5 kg, total weight 50 kg
- All pulleys are fitted with ball bearings
- The rope length can be adjusted via a serial perforation irrespective of the rope height
- Frame: silver-anthracite

Item no. 48125 W x D x H: 31 x 34 x 232 cm

PINOFIT® Vertical pulley with latissimus bracket

- Load: 5 x 2 kg, 10 x 5 kg, total weight 60 kg
- Direct power transmission
- The rope length can be adjusted via a serial perforation
- Frame: silver-anthracite

Item no. 48134 W x D x H: 31 x 89 x 236 cm

Latissimus bracket for pulley apparatus (not shown)

Item no. 48128 Aluminium, 104 cm wide

Tripod for pulley apparatus or vertical pulley

Item no. 48129

WALL BARS PRO

The multifunctional training station.

MADE
IN
GERMANY



Pictured optional equipment available at extra cost: metal ceiling grille and adjustable stainless steel tubes.

Shipping costs see p. 176

PINOFIT® wall bars Pro

- A must for every practice
- With manifold extension options
- Extremely durable metal side boards
- 11 round rungs made of beech plywood according to DIN regulations, incl. 1 overhang rung
- Distance between rungs approx. 17.5 cm
- Incl. 6 wall brackets
- Supplied pre-assembled
- Approx. 240 x 100 cm

Item no. 41640



✓ MEDICAL PRODUCT CE



Dip bar made from aluminium

- To be inserted in the bars – with arm and back cushion in dark grey

Item no. 41646
Approx. W x D: 71 x 62 cm



Roof grid from metal

- ideal for sling therapy & extensions
- Uncomplicated to screw on
- Incl. mounting accessories

Item no. 41647
Approx. W x D: 90 x 80 cm



Stainless steel pull-up bar

- To be inserted in the wall bar
- Versatile applications

Item no. 41648
Approx. W x D: 104 x 50 cm



Two stainless steel adjustment tubes

- Each with two ports for attaching gymnastic bands, for example
- Ideal for a variety of exercises
- Bars with numbering

Item no. 41644
Approx. 240 cm height



Supplied without pull-up bar.



Wall mount for accessories Wall bars Pro

- For inserting pull-up bar and dip bar (see p. 118)
- For practical storage of the accessories for the wall bars PRO

Item no. 41651 W x D x H: 52 x 9 x 20 cm



MADE
IN
GERMANY

Wall bars

- End piece from robust finger-jointed aspen wood
- Rungs from beech plywood with round profile
- Distance between rungs approx. 15 cm
- Rungs according to DIN regulations
- 4 wall brackets
- Supplied assembled
- Available in 2 sizes
- Including 1 overhang rung
- Including brackets for brick walls
- Made in Germany

Item no. 41576 Special physio version, approx. 230 x 100 cm, 12 round rungs

Item no. 41575 Approx. 210 x 80 cm, 11 round rungs



High quality slings

- Durable cotton quality
- Extra stable stitching
- Colour: cotton nature
- Reinforced stainless steel ring suspension
- Perfect for sling therapy



i High quality slings in various designs can be found online.

MADE
IN
GERMANY

COMPLETE PHYSIO PACKAGE

To improve strength, flexibility & coordination.

MADE
IN
GERMANY



✓ MEDICAL PRODUCT CE

- 1 PINOFIT® functional presses
- 2 PINOFIT® pulley apparatus standard wall models including hand grips
- 1 hand cuff (not shown)
- 1 foot cuff (not shown)
- 1 PINOFIT® vertical pulley wall model including latissimus bracket
- 1 PINOFIT® training bench
- 1 PINOFIT® angle table

Item no. 48149

i Detailed information and all articles also available individually online.

TRAINING ROOM ACCESSORIES



Swedish apparatus

- Pine wood treated with beeswax
- Smooth with balance beam and mounting fixture for the wall bar

Item no. 41377
Approx. length 200 cm x height 40 cm

Trampoline Trimlin-Med Plus

- Maximum load 120 kg
- TÜV/GS-tested
- Foldable feet with floor protectors
- High quality polypropylene trampoline with safety cover

Item no. 44511
ø approx. 102 cm, height approx. 25 cm

PINO correction mirror

- Mirror surface made from crystal glass
Approx. 58 x 149 cm for a clear and brilliant image
- Frame made from powder-coated sheet steel
- Can be tilted backwards and forwards
- Low rolling resistance wheels for mobile use
- Adjustable rollers in front

Item no. 48282
W x D x H: 60 x 70 x 166 cm



PINOFIT® mat trolley

- For up to 20 mats 60 cm wide
- Powder-coated sheet steel
- Extremely robust and durable
- Low rolling resistance wheels for mobile use
- Supplied not assembled

Item no. 48281 W x D x H: 132 x 70 x 107 cm



PINO wall mount fixture

- For 5-10 mats, supplied without fixtures and mats
- Silver powder-coated sheet steel

Item no. 48194 W x D x H: 135 x 25 x 20 cm




PINO wall mount fixture for balls

- For three balls
- Supplied without fixtures or balls
- Silver powder-coated sheet steel

Item no. 48193 W x D x H: 198 x 45 x 20 cm

THERAPY BEDS & PRACTICE STOOLS





TOP QUALITY MADE IN GERMANY



PINO therapy beds are made **exclusively in Germany**. The entire development process takes place in-house, from the product idea to the finished production model. The products are customized in **our own production** facilities with loving attention to detail and the highest levels of quality, workmanship and durability. All in all, reliability is the top priority. With our years of experience in the field of physiotherapy, we meet the highest standards of quality, workmanship and durability. To meet your needs and requirements, we work hand in hand with experienced physiotherapists.

✓ MEDICAL PRODUCT CE

Which bed is right for me?

				
Technical details	PINOTREND	PINOPRO	PINO XPERT+	PINO Q
	The entry-level model	The all-rounder with extra configuration options	Narrow frame for ergonomically optimised work on the patient	Narrow design without frame
Lying surface	Standard	Standard	Shape-optimised	Shape-optimised
Available bed widths	65, 70, 80 cm	65, 70, 80 cm	65, 80 cm	65, 80 cm
Available lying surfaces		2-part lying surface with adjustable footrest*	2-part lying surface with adjustable footrest*	2-part lying surface with adjustable footrest*
		2-part lying surface with electrical flexion**	2-part lying surface with adjustable headrest**	
		Thermal therapy bed	Thermal therapy bed	
Available height adjustments	Foot switch	Foot switch	Foot switch	Foot switch
	Hand lever	All-round control bar Hand lever	All-round control bar Hand lever	Height adjustment in frame, hand lever
Position can be shifted	Wheel axis lifting system	Central wheel lifting system	Integrated central wheel lifting system as standard	Two wheels integrated in the frame as standard (one axle)
Cushioned	Single-stitched	Single-stitched	Double-stitched	Single-stitched
Load capacity	200 kg	250 kg	250 kg	250 kg

* The 2-piece bed surface with hinged foot section enables the perfect elevated positioning

**The flexion position allows shock and extension treatments for ideal spinal positioning



INNOVATIVE HEADRESTS

The right solution for every treatment.



Standard Single piece headrest

The headrest provides all the basic functions for performing professional physiotherapy treatments. Thanks to the oval nose slot, the model is also suitable for treatments in the prone position. You can adjust the inclination continuously from approx. -90° to +30°.



Four-piece headrest Move

In addition to the supine and prone positions, the Move model is used for the fixation of extremities in the treatment of the cervical spine, lower legs, knees, arms and hands. With the help of the sliding face opening, you can ideally adjust the distance to suit your patient's specified head size. Additional armrests allow for an even more comfortable ergonomic positioning of patients. The continuously adjustable inclination ranges from approx. -60° to +35°.



Extra-long headrest with adjustable backrest (only for therapy bed XPERT+)

Thanks to the integrated nose slot, this long headrest allows you to lie on your stomach as well as on your back. A raised-knees position is also possible. For comfortable treatment while sitting, adjust the back section continuously up to +60° - simply by using the lever.



TOP SELLER

Three-piece headrest Flex

In addition to the supine and prone positions, the Flex model is also used for treatments in the sitting position. A raised-knees position is also possible. Flexibility is the focus here: adjust the height of both the headrest and armrests and tilt the headrest from -90° to +35°.

Three-piece headrest

Just like the 1-piece headrest, the 3-piece headrest is recommended for the supine and prone position. In addition, it offers height-adjustable and swivelling armrests. The tilt angle of the headrest ranges from approx. -90° to +35° (-45° to +35° for PINO Xpert+ and PINO Q) and can also be adjusted continuously and smoothly via a gas pressure spring.



TOP SELLER

Three-piece headrest FIX

Just like the 1-piece headrest, the Fix headrest is recommended for the supine and prone position. The Fix version, however, additionally offers tilt- and height-adjustable armrests. The tilt angle of the headrest ranges from approx. -55° to +35° and can also be adjusted continuously and smoothly via a gas pressure spring.



Our highlights for you

- ⊕ Arrange an appointment for a visit by phone with your PINO customer consultant, our hotline or consultation via video conference: www.pinoshop.de/digitale-beratung
- ⊕ Let us advise you - visit our showroom in Büren

HANDMADE – PRODUCED IN-HOUSE

The highest quality guaranteed

- Very stable and continuous cross frames with **5-year guarantee**
- Exclusive use of German quality **motors with 2-year guarantee**
- The high-quality upholstery is of above-average strength and consists of a combination of various foam materials with different degrees of hardness for optimum comfort and ultimate shape retention
- Extensive range of upholstery and frame colours for unlimited opportunities to create your individual choice of therapy bed

Quality assurance

- Regular **tests of long-term usage** - prototypes of each new model are subjected to a minimum of **22,000 full lifts** (complete up and down movements) with a continuous load of **160 kg**
- Motors, are hand-made and subjected to a maximum load test after production.



PINO XPERT+


Sophisticated design and maximum functionality.




Pictured optional equipment available at extra cost: all-round control bar for practical height adjustment.



MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT | SUPPORTS A WEIGHT OF UP TO 250 KG | MEDICAL PRODUCT CE



Two-part lying surface with adjustable footrest



Single-part, extra-long headrest

PINO XPERT+ therapy bed

- Standard control with electric height adjustment from 52 to 96 cm in just 16 sec.
- Control via low-maintenance German quality motor
- Two-tone frame in silver and dark grey, tidy look with hidden cables
- Ergonomic work on patient due to tapered foot area, load capacity 250 kg
- With hand-sewn premium upholstery 7 cm thick with antimicrobial synthetic leather cover
- Shape-optimised, one-piece headrest with nose slot, continuously and smoothly adjustable from -90° to +30°
- The robust misalignment-proof cross frame is electrical and is fitted with maintenance-free dry plain bearings
- **With 5-year warranty for the frame and 2-year warranty for the motor**
- The safety precautions of the PINO therapy beds comply with the current specifications of the BfArM (German Federal Institute for Drugs and Medical Devices)

Item no. 48080



Standard highlights:

- ⊕ Double-stitched premium padding
- ⊕ Central wheel lifting system
- ⊕ No visible cables
- ⊕ Ergonomically optimised lying surface
- ⊕ German quality motor
- ⊕ Frame in two-tone dark/silver grey
- ⊕ Height adjustable in only 16 sec.

Select your optional extras here:
(continued on p. 143)



Foot-operated lever instead of hand-operated lever

Item no. 48154




Oval nose slit

Item no. 48121



Comfort padding height 10 cm

Item no. 48454



Padding width 80 cm

Item no. 48453



Four-part headrest MOVE
65 cm wide, approx. -60°/+35° continuously adjustable via gas pressure spring

Item no. 48386



Three-part headrest
65 cm wide, -45°/+35° continuously adjustable via gas pressure spring, with additional, height-adjustable armrests

Item no. 48383



Three-part headrest FIX
With tiltable armrests, 65 cm wide, freely adjustable via gas pressure spring

Item no. 48384



Three-part headrest FLEX
65 cm wide, -90°/+35° continuously adjustable via gas pressure spring, perfect for treatments in supine and sitting positions, with additional, height-adjustable armrests

Item no. 48385



Single-part, extra-long headrest With hinged back
65 cm wide, 75 cm long, +60° continuously adjustable via gas pressure spring, with integrated oval nose slit for treatments in the prone position

Item no. 48388



Select your optional extras here:




Two-part lying surface with adjustable footrest
 Leg surface +75° is conveniently continuously adjustable via gas pressure spring

Item no. 48380




Two-part lying surface with adjustable backrest
 Headrests +60° conveniently continuously adjustable via gas pressure spring

Item no. 48381



Heated lying surface
 Bed loses status as medical device, only available for 1-piece lying surface

Item no. 48382



ALL-ROUND CONTROL BAR
 for practical height adjustment of the therapy bed from any position

Item no. 48452



Protective mat
 Approx. 60 x 65 cm, made from durable and wipeable PVC material

Item no. 40845

Pictured optional equipment available at extra cost: all-round control bar for practical height adjustment.

PINOTREND


The powerful starter model.



MADE OF ANTIMICROBIAL IMITATION LEATHER
DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF
UP TO 200 KG

MEDICAL PRODUCT CE



Standard highlights:

- ⊕ German quality motor
- ⊕ 7 cm comfort padding

PINOTREND therapy bed

- With standard electrical height adjustment via a manual lever, from 48 to 96 cm in only 25 seconds,
- 1-part surface 200 x 65 cm, with maximum load of 200 kg
- Headrest with nose slit – 90°/+30° freely adjustable via a gas pressure spring
- Extra stable misalignment-proof cross frame
- Equipped with maintenance-free dry sliding bearings
- 7 cm comfort upholstery with antimicrobial synthetic leather cover
- Rounded, reinforced corners, without seams
- **5-year guarantee for the frame and 2-year guarantee for the motor**
- The safety precautions of the PINO therapy beds comply with the current specifications of the BfArM (German Federal Institute for Drugs and Medical Devices)

Item no. 48100

Item no. 48101 Hydraulic with foot lever

Select your optional extras here:

<p>Turbo motor for height adjustment in only 16 seconds <i>not available for PINOPRO hydraulic</i> Item no. 49150</p>	<p>Foot-operated lever instead of hand-operated lever <i>Not available for PINOPRO hydraulic</i> Item no. 48154</p>
<p>Padding width 70 cm Item no. 48415</p>	<p>Padding width 80 cm Item no. 48413</p>
<p>Comfort padding height 10 cm Item no. 48189</p>	<p>Oval nose slit Item no. 48121</p>

Select your optional extras here:



Four-part headrest MOVE
65 cm wide, approx. -60°/+35° continuously adjustable via gas pressure spring

Item no. 48301



Three-part head piece
65 cm wide, -90°/+35° continuously adjustable via gas pressure spring, with additional, height-adjustable armrest

Item no. 48300



Three-part headrest FIX
With tiltable armrests, 65 cm wide, freely adjustable via gas pressure spring

Item no. 48307



Three-part headrest FLEX
65 cm wide, -90°/+35° continuously adjustable via gas pressure spring, perfect for treatments in supine and sitting positions, with additional, height-adjustable armrests

Item no. 48327



Wheel axle lifting system
Comfortable starter mechanics via two easy moving foot levers for the control of the wheels

Item no. 48153




Paper roll holder
Holder integrated into the frame at the foot end

Item no. 48152



Paper rolls
9 rolls of 50 m each, 2-ply, perforated, 143 sheets, 59 cm wide

Item no. 40890



Protective mat
Approx. 60 x 65 cm, made from durable and wipeable PVC material

Item no. 40845

PINOPRO

The all-rounder for the highest requirements.






Pictured optional equipment available at extra cost: 2-part lumbar section incl. electrical flexion, central wheel-lifting system, four-part headrest and all-round control bar.



MADE OF ANTIMICROBIAL IMITATION LEATHER
DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 250 KG

MEDICAL PRODUCT CE



PINOPRO therapy bed

- With standard electrical height adjustment via a manual lever, from 48 to 96 cm in only 25 sec.
- 1-part surface 200 x 65 cm, with maximum load of 250 kg
- Headrest with nose slit $-90^{\circ}/+30^{\circ}$ freely adjustable via a gas pressure spring
- Extra robust misalignment-proof frame for maximum stability ensured by large round tubes
- Equipped with maintenance-free dry sliding bearings
- 7 cm comfort padding with antimicrobial synthetic leather cover
- Rounded, reinforced corners, without seams
- **5-year guarantee for the frame and 2-year guarantee for the motor**
- The safety precautions of the PINO therapy beds comply with the current specifications of the BfArM (German Federal Institute for Drugs and Medical Devices)

Item no. 48110


Item no. 48111 Hydraulic with foot lever



Standard highlights:

- ⊕ German quality motor
- ⊕ 7 cm comfort padding
- ⊕ Particularly sturdy frame can withstand loads of up to 250 kg
- ⊕ Wide range of additional equipment for all types of therapy

Select your optional extras here:
(continued on p. 148)



Turbo
motor for height adjustment in only 16 seconds
not available for PINOPRO hydraulic

Item no. 49150




Foot-operated lever instead of hand-operated lever
Not available for PINOPRO hydraulic

Item no. 48154




Comfort padding height 10 cm

Item no. 48189



Padding width 70 cm

Item no. 48415



Padding width 80 cm

Item no. 48413



Oval nose slit

Item no. 48121



Four-part headrest MOVE
65 cm wide, approx. $-60^{\circ}/+35^{\circ}$ continuously adjustable via gas pressure spring

Item no. 48301



Three-part head piece
65 cm wide, $-90^{\circ}/+35^{\circ}$ continuously adjustable gas pressure spring, with additional, height-adjustable armrests

Item no. 48300



Three-part headrest FIX
With tiltable armrests, 65 cm wide, freely adjustable via gas pressure spring

Item no. 48307

Select your optional extras here:




Three-part headrest FLEX
65 cm wide, -45°/+35° continuously adjustable via gas pressure spring, perfect for treatments in supine and sitting positions, with additional, height-adjustable armrests
Item no. 48327



Two-part lying surface
With safety distance between bed padding and frame, leg surface +75° comfortably infinitely adjustable via gas pressure spring
Item no. 48305



Two-part lying surface with electrical flexion
Not available for PINOPRO hydraulic
Item no. 48310



Thermal therapy bed
Bed loses status as medical product
Not available for PINOPRO hydraulic
Item no. 48104



Central wheel lifting system
Especially comfortable starter mechanics via just one easy moving foot lever for the control of all four wheels
Item no. 48311



All-round control bar
For practical height adjustment of the therapy bed from any position
Not available for PINOPRO hydraulic
Item no. 48312



Paper roll holder
Holder integrated into the frame at the foot end
Item no. 48152



Paper rolls
9 rolls of 50 m each, 2-ply, perforated, 143 sheets, 59 cm wide
Item no. 40890



Protective mat
Approx. 60 x 65 cm, made from durable and wipeable PVC material
Item no. 40845

PINO BASIC

The reliable basic model.



MADE OF ANTIMICROBIAL IMITATION LEATHER
DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 200 KG

MEDICAL PRODUCT CE

PINO BASIC therapy bed

- 75 cm high
- 2-part surface 200 x 65 cm, with maximum load of 200 kg
- Headrest with nose slit -90°/+30° freely adjustable via a gas pressure spring
- 7 cm comfort padding with antimicrobial synthetic leather cover
- Rounded, reinforced corners, without seams
- 5-year guarantee for the frame

Item no. 48130

Padding width 70 cm
Item no. 48415

Comfort padding height 10 cm
Item no. 48189

Padding width 80 cm
Item no. 48413

Oval nose slit
Item no. 48121

COMFORT PADDING FOR TREATMENT BEDS

OFFER YOUR CUSTOMERS FIRST-CLASS ADDITIONAL EXTRA COMFORT

Comfort padding for treatment beds

- Extra soft lying surface with super soft polyurethane cover
- Hand-sewn luxury upholstery 9 cm thick
- Upholstery cover has rounded reinforced corners and is water-proof, resistant to oils and greases and provides protection from bacteria
- Heated lying surface for an additional charge, for even more treatment comfort
- Rubberised underside prevents the support from slipping on the lying surface

Comfort padding for treatment beds
Item no. 40747

Comfort padding for treatment beds, thermal
Item no. 40748



PINO Q – INTELLIGENTLY SIMPLE

With a clever design and perfect functionality.



Optional extras for an additional charge: 2-part lumbar area and 3-part headrest.





ENTSPRICHT
DEN AKTUELLEN
SICHERHEITSEMPFEHLUNGEN
DES BfArM

MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 250 KG

MEDICAL PRODUCT CE



PINO Q
Two-part lumbar section, item no. 48106

PINO Q THERAPY BED

- Standard equipment: electrical height adjustment via a frame integrated foot lever between 48 and 96 cm (in only 16 seconds) and a two-part surface
- Thanks to the wheels incorporated in the front frame, the PINO Q can be easily moved to a different location
- The head piece with breathing hole can be adjusted continuously by $-45^{\circ}/+30^{\circ}$ via a gas pressure spring
- Shape-optimised 1-part lying surface 200 x 65 cm for easier work on the patient,
- max. load capacity 250 kg
- Rounded, reinforced corners, without seams
- Stable, misalignment-proof cross frame with maintenance-free dry plain bearings
- **With 5-year-warranty for the frame and 2-year-warranty for the motor**
- The safety precautions of the PINO therapy beds comply with the current specifications of the BfArM (German Federal Institute for Drugs and Medical Devices)

Item no. 48109



Standard highlights:

- ⊕ German quality turbo motor
- ⊕ Ergonomically optimised lying surface
- ⊕ For a weight up to 250 kg
- ⊕ No visible cables

Select your optional extras here:



Four-part headrest MOVE
65 cm wide, approx. $-60^{\circ}/+35^{\circ}$ continuously adjustable
via gas pressure spring

Item no. 48323



Three-part headrest
65 cm wide, $-45^{\circ}/+35^{\circ}$ continuously adjustable
via gas pressure spring
, with additional, height-adjustable armrests

Item no. 48321



Three-part headrest FIX
With tiltable armrests,
65 cm wide, $-55^{\circ}/+35^{\circ}$ freely adjustable via gas
pressure spring

Item no. 48333



Three-part headrest FLEX
65 cm wide, $-45^{\circ}/+35^{\circ}$ continuously adjustable
via gas pressure spring, perfect for treatments
in supine and sitting positions, with additional,
height-adjustable armrests

Item no. 48334



Two-part lying surface
Leg area $+75^{\circ}$ freely adjustable via gas pressure
spring

Item no. 48106



Protective mat
Approx. 60 x 65 cm, made from durable
and wipeable PVC material

Item no. 40845

Additional hand-operated lever

Item no. 48119

Additional height adjustment within the frame

Item no. 48292

Comfort padding height 10 cm

Item no. 48189

Additional foot-operated lever

Item no. 48298

Padding width 80 cm

Item no. 48416

Oval nose slit

Item no. 48121

PINOSPA®

Exclusive design for the highest standards.



Pictured optional equipment available at extra cost: all-round control bar, central wheel-lifting system.




ENTSPRICHT DEN AKTUELLEN SICHERHEITSEMPFEHLUNGEN DES BfArM

MADE OF ANTIMICROBIAL IMITATION LEATHER
DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 200 KG

MEDICAL PRODUCT CE



Plastic surface in oak look optional



Mocha



Nature



White

PINOSPA® therapy bed

- Height adjustable as standard with new-style lifting column turbo motor in only 16 seconds, including lockout box
- Very quiet and smooth
- Incl. magnetic hand switch
- 1-part surface 200 x 70 cm, with maximum load of 200 kg
- 1-part headrest with nose slit $-90^{\circ}/+30^{\circ}$ freely adjustable via a gas pressure spring
- 7 cm comfort upholstery with antimicrobial synthetic leather cover
- Durable plastic surface in exclusive oak look, easy to disinfect
- Rounded, reinforced corners, without seams
- **2-year-warranty on the motor**
- The safety precautions of the PINO therapy beds comply with the current specifications of the BfArM (German Federal Institute for Drugs and Medical Devices)

Item no. 48095



Standard highlights:

- ⊕ German quality turbo motor
- ⊕ Frame in oak look
- ⊕ No visible cables
- ⊕ 70 cm padding width

Select your optional extras here:



Four-part headrest MOVE
65 cm wide, approx. $-60^{\circ}/+35^{\circ}$ continuously adjustable via gas pressure spring

Item no. 48301



Three-part headrest
65 cm wide, $-90^{\circ}/+35^{\circ}$ continuously adjustable gas pressure spring, with additional, height-adjustable armrests

Item no. 48300



Three-part headrest FLEX
65 cm wide, $-45^{\circ}/+35^{\circ}$ continuously adjustable via gas pressure spring, perfect for treatments in supine and sitting positions, with additional, height-adjustable armrests

Item no. 48327



Central wheel lifting system
Especially comfortable starter mechanics via an easy moving foot lever for the control of all four wheels

Item no. 48093



All-round control bar
For practical height adjustment of the therapy bed from any position

Item no. 48094



Protective mat
Approx. 60 x 65 cm, made from durable and wipeable PVC material

Item no. 40845

Thermal therapy bed
Bed loses status as medical device

Item no. 40105

Height adjustment
from 52 to 72 cm

Item no. 44051

Height adjustment
from 62 to 82 cm

Item no. 44052

Foot-operated lever instead of hand-operated lever

Item no. 48154

Oval nose slit

Item no. 48121

Padding width 80 cm

Item no. 48413

Comfort padding height 10 cm

Item no. 48189

PINOSPA® SOFT

New standards for spa and wellness treatments.



MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 200 KG

MEDICAL PRODUCT CE

PINOSPA® Soft therapy bed

- Extra soft lying surface with super soft polyurethane cover
- Hand-sewn luxury upholstery 10 cm thick
- Upholstery cover has rounded reinforced corners and is water-proof, resistant to oils and greases and provides protection from bacteria
- Extra-long headrest with nose slit and adjustable back, infinitely adjustable from 0° to +80° by means of a gas pressure spring
- Electrically height adjustable via magnetic hand switch in just 16 seconds
- Two different height adjustments freely selectable: 52–72 cm and 62–82 cm
- Frame with durable plastic surface, particularly easy to disinfect
- L x W x H: 200 x 82 x 52–72 cm or 62–82 cm
- **2-year-warranty on the motor**
- The safety precautions of the PINO therapy beds comply with the current specifications of the BfArM (German Federal Institute for Drugs and Medical Devices)

Item no. 48086

Select your optional extras here:

Foot-operated lever instead of hand-operated lever
Item no. 48154

All-round control bar
Item no. 48094

Central wheel lifting system
Item no. 48093

Oval nose slit
Item no. 48121

Thermal therapy bed
Bed loses status as medical device
Item no. 40105

PINOPRO BOBATH


The specialist model for Bobath therapy.



MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 250 KG

MEDICAL PRODUCT CE



PINOPRO Bobath special fittings: single piece headrest without nose slit item no. 48161

PINOPRO Bobath therapy bed

- With standard electrical height adjustment via a manual lever, from 48 to 96 cm in only 25 sec.
- Single-part surface 200 x 100 cm, max. weight 250 kg
- Extra stable misalignment-proof frame for maximum stability ensured by large round tubes
- Equipped with maintenance-free dry sliding bearings
- 7 cm comfort padding with antimicrobial synthetic leather cover
- Rounded, reinforced corners, without seams
- **5-year guarantee for the frame and 2-year guarantee for the motor**
- The safety precautions of the PINO therapy beds comply with the current specifications of the BfArM (German Federal Institute for Drugs and Medical Devices)

Item no. 48115

Select your optional extras here:

Turbo motor
For height adjustment in only 16 seconds
Item no. 49150

1-piece headrest length 34 cm
*Without nose slit
Lying surface reduces to 161 cm*
Item no. 48161

Comfort padding height 10 cm
Item no. 48190

Padding width 120 cm
Item no. 48147

Foot-operated lever instead of hand-operated lever
Item no. 48154

Central wheel lifting system
Item no. 48088

All-round control bar
Item no. 48312

PINOPRO LYMPH

The special model for lymph treatments.





MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT | SUPPORTS A WEIGHT OF UP TO 250 KG | MEDICAL PRODUCT CE

PINOPRO LYMPH therapy bed

- With standard electrical height adjustment via a manual lever, from 48 to 96 cm in only 25 sec.
- 5-part surface 200 x 65 cm, with maximum load of 250 kg
- 4-part foot part whereby each part can be individually and infinitely adjusted via a lockable gas pressure spring
- Headrest with nose slit $-90^{\circ}/+30^{\circ}$ freely adjustable via a gas pressure spring
- Extra robust misalignment-proof frame for maximum stability ensured by large round tubes
- Equipped with maintenance-free dry sliding bearings
- 7 cm comfort padding with antimicrobial synthetic leather cover
- Rounded, reinforced corners, without seams
- **5-year guarantee for the frame and 2-year guarantee for the motor**
- The safety precautions of the PINO therapy beds comply with the current specifications of the BfArM (German Federal Institute for Drugs and Medical Devices)

Item no. 48135 Electric with manual switch

Item no. 48132 Hydraulic with foot lever



Select your optional extras here:



Four-part headrest MOVE
65 cm wide, approx. $-60^{\circ}/+35^{\circ}$ continuously adjustable via gas pressure spring

Item no. 48301



Three-part headrest
65 cm wide, $-90^{\circ}/+35^{\circ}$ continuously adjustable gas pressure spring, with additional, height-adjustable armrests

Item no. 48300



Three-part headrest FIX
With tiltable armrests, 65 cm wide, $-55^{\circ}/+35^{\circ}$ freely adjustable via gas pressure spring

Item no. 48307



Three-part headrest FLEX
65 cm wide, $-45^{\circ}/+35^{\circ}$ continuously adjustable via gas pressure spring, perfect for treatments in supine and sitting positions, with additional, height-adjustable armrests

Item no. 48327



Central wheel lifting system
Especially comfortable starter mechanics via just one easy moving foot lever for the control of all four wheels

Item no. 48311



All-round control bar
For practical height adjustment of the therapy bed from any position
Not available for Pinopro Lymph hydraulic

Item no. 48312



Paper roll holder
Holder integrated into the frame at the foot end

Item no. 48152

Paper rolls
9 rolls of 50 m each, 2-ply, perforated, 143 sheets, 59 cm wide

Item no. 40890

Turbo motor
For height adjustment in only 16 seconds
Not available for Pinopro Lymph hydraulic

Item no. 49150

Foot-operated lever instead of hand-operated lever
Not available for Pinopro Lymph hydraulic

Item no. 48154

Oval nose slit

Item no. 48121

Padding width 80 cm

Item no. 48418

Comfort padding height 10 cm

Item no. 48189

PRACTICE STOOL COMPACT

With round seat and hand ring.

New



Particularly narrow base, for working close to the patient.

PRACTICE STOOL

With round seat and hand ring.



MADE IN GERMANY

PRACTICE STOOL

With round seat and extra wide armrest and backrest.

New



MADE IN GERMANY



Extra-wide arm and backrest: resting forearms allows calmer movements during treatment.

5 5-YEAR-WARRANTY FOR THE FRAME

MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 100 KG

5 5-YEAR-WARRANTY FOR THE FRAME

MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 100 KG

- Round, padded seat
- Comfortable height adjustment via ring underneath the seat
- Comfort upholstery with antimicrobial synthetic leather cover
- With base made of black, robust plastic as standard
- Standard casters for hard flooring
- Upholstery colours are free to choose; please see cover

Item no. 48376 Seat ø 36 cm, height adjustable from 52–71 cm


- Round, padded seat
- Comfortable height adjustment via ring underneath the seat
- Comfort upholstery with antimicrobial synthetic leather cover
- Cross-shaped base in black powder-coated aluminium as standard
- Standard casters for hard flooring
- Upholstery colours are free to choose; please see cover

Item no. 48215 Seat ø 36 cm, height adjustable from 52–71 cm

- Round, padded seat
- Comfortable height adjustment via ring underneath the seat
- Armrest continuously height adjustable 78–85 cm
- Comfort upholstery with antimicrobial synthetic leather cover
- Cross-shaped base in black powder-coated aluminium as standard
- Standard casters for hard flooring
- Upholstery colours are free to choose; please see cover

Item no. 48375 Seat ø 36 cm, height adjustable from 52 to 71 cm

Choose your favourite accessories for your practice stool



Dynamic 360° seat mechanics
with hand-operated lever

Item no. 48296



Cross-shaped base
in polished aluminium
Only for stool 48215

Item no. 48288



Caster accessories
Wheels are blocked under weight

Item no. 48287

Choose your favourite accessories for your practice stool



Cross-shaped base
in polished aluminium

Item no. 48288



Caster accessories
Wheels are blocked under weight

Item no. 48287

PRACTICE STOOL

With saddle seat and hand ring.



PRACTICE STOOL

With saddle seat and hand lever.



MADE IN GERMANY

5 5-YEAR-WARRANTY FOR THE FRAME

MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 100 KG

- Padded saddle seat
- Comfortable height adjustment via ring underneath the seat
- Comfort upholstery with antimicrobial synthetic leather cover
- Cross-shaped base in black powder-coated aluminium as standard
- Standard casters for hard flooring
- Upholstery colours are free to choose; please see cover

Item no. 48230 Seat ø 36 cm, height-adjustable from 52 to 71 cm

Choose your favourite accessories for your practice stool



Caster accessories
Wheels are blocked under weight

Item no. 48287



Cross-shaped base
in polished aluminium

Item no. 48288



Dynamic 360° seat mechanics
with hand-operated lever
Only for saddle seat stool!

Item no. 48296

- Padded seat in saddle form
- Comfortable height adjustment via hand lever
- Comfort upholstery with antimicrobial synthetic leather cover
- Cross-shaped base in black powder-coated aluminium as standard
- Standard casters for hard flooring
- Upholstery colours are free to choose; please see cover

Item no. 48220 Seat ø 36 cm, height adjustable from 52 to 71 cm

PRACTICE STOOL

With round seat and backrest.

- With extra wide, height-adjustable backrest for more comfort and an upright and stable sitting position
- Comfort upholstery with antimicrobial synthetic leather cover
- Cross-shaped base in black powder-coated aluminium as standard
- Casters for hard flooring
- Upholstery colours are free to choose; please see cover


Item no. 48240 Seat ø 36 cm, height-adjustable from 52 to 71 cm

Choose your favourite accessories for your practice stool



Cross-shaped base
in polished aluminium

Item no. 48288



Caster accessories
Wheels are blocked under weight

Item no. 48287



MADE IN GERMANY

5 5-YEAR-WARRANTY FOR THE FRAME

MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 100 KG



MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 100 KG SUPPORTS A WEIGHT OF UP TO 100 KG SUPPORTS A WEIGHT OF UP TO 200 KG

Stacking stool

- Particularly robust and versatile
- Free choice of upholstery colour
- Maximum load 100 kg
- Material: birch wood

Item no. 48270
Seat ø 37 cm, W x H 45 x 46 cm

Stool with four legs

- made from rectangular tubing
- Comfort upholstery with antimicrobial synthetic leather cover
- Colour of upholstery freely selectable

Item no. 48227
Size approx. 40 x 40 cm, height 47 cm

Stool with four legs PRO

- made from rectangular tubing
- Comfort upholstery with antimicrobial synthetic leather cover
- Colour of upholstery freely selectable

Item no. 48365
Size approx. 40 x 40 cm, height 47 cm

PINO MASSAGE BENCH MOBIL

Easy handling for mobile therapy.



Headrest with arm sling for a more comfortable prone position.
Height individually adjustable.



MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 180 KG



PINO MOBIL Profi portable massage bench alu

Including carrier bag and arm sling

- Surface approx. 183 x 62 cm
- Frame made from extra light aluminium
- Headrest can be adjusted positively and negatively, with arm sling
- Height adjustment from approx. 60 to 81 cm
- Suitable for a weight of up to 180 kg, tare weight approx. 12 kg
- Extra strong padding, 10 cm, for very high comfort
- Including carrier bag and arm sling

Mint

Item no. 42283 Alu. frame

Dark grey

Item no. 42288 Alu. frame

Sky blue

Item no. 42281 Alu. frame

Lime

Item no. 42287 Alu. frame

Beige

Item no. 42282 Alu. frame

PINO MOBIL Profi portable massage bench

Including carrier bag, arm sling and semi roll

- Surface approx. 183 x 76 cm
- Frame from high quality beech, varnished
- Headrest can be adjusted positively and negatively, with arm sling
- Height adjustment from approx. 60 to 86 cm
- Suitable for a weight of up to 180 kg, tare weight circa 18 kg
- Extra strong padding, 11 cm, for very high comfort
- Including carrier bag, arm sling and semi roll

Mint

Item no. 42253 Wood frame

Dark grey

Item no. 42258 Wood frame

Sky blue

Item no. 42254 Wood frame

Lime

Item no. 42257 Wood frame

Beige

Item no. 42256 Wood frame



MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

SUPPORTS A WEIGHT OF UP TO 250 KG



PINO MOBIL Profi portable massage bench oval

Including carrier bag and arm sling

- Simple folding function with high quality closures
- High lying comfort thanks to soft foam padding
- Easy-care imitation leather cover with strong stitching
- Padding thickness: 6 cm
- Two-part lying surface approx. 195 x 76 cm
- Back section can be set up to 70° in increments
- Headrest can be inserted at the front and back of the bench
- Headrest and head padding can be stowed inside the couch during transport
- Sturdy beech wood frame, height adjustable via two screws from approx. 61 to 81.5 cm
- Suitable for a weight of up to 250 kg

Dark grey

Item no. 42276 Wood frame

Lime

Item no. 42271 Wood frame

PINO MOBIL Profi portable massage bench alu extra-wide

Including carrier bag and arm sling

- Simple folding function
- High lying comfort thanks to soft foam padding
- Easy-care imitation leather cover with strong stitching
- Padding thickness: 6 cm
- Two-part lying surface: 184 x 71 cm
- Back/foot section can be set up to 70° in increments
- Headrest can be inserted at the front and back of the bench
- Headrest and head padding can be stowed inside the couch during transport
- Sturdy aluminium frame, height adjustable from 61 to 81.5 cm
- Suitable for a weight of up to 250 kg

Dark grey

Item no. 42274 Alu. frame

Lime Alu. frame

Item no. 42273 Alugestell

PINO MASSAGE CHAIR PRO

New version of our classic.



PINO MOBIL Profi massage stool

Incl. rolling case and chest pad

- Ergonomic design
- Comfortable upholstery with easy-care imitation leather cover
- Headrest, armrest, chest and seat padding individually adjustable
- Headrest adjustable in height and inclination, with removable U-cushion
- Optimal patient positioning for any body size
- Stable and secure up to a weight of 200 kg
- Easy transport in rolling case

Dark grey

Lime

Item no. 48369

Item no. 48339

 MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

 SUPPORTS A WEIGHT OF UP TO 200 KG



Chest support



Adjustable front frame



Adjustable headrest



Adjustable chest pad



PINO MASSAGE CHAIR MOBIL

The mobile massage support.



PINO MOBIL Profi massage stool

Including carrier bag and chest padding

- Extremely light aluminium frame
- Individually adjustable for differing heights, headrest with removable U-pillow
- Height and incline of the headrest are adjustable
- Stable and secure up to a weight of 150 kg
- Extremely light, tare weight 8.8 kg
- Padding from easy-care artificial leather
- Including carrier bag and removable chest padding

Grey

Lime

Item no. 48338

Item no. 48337

 MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

 SUPPORTS A WEIGHT OF UP TO 150 KG



For relaxed standing during work

PINO anti-fatigue mat

- Ergonomic mat for work standing, e.g. work at the therapy bed
- The optimum cushioning properties and special surface structure are designed to reduce signs of fatigue after prolonged work standing up
- Anti-slip with nubbed surface, single-layer, durable and shock absorbing
- Suitable for dry work areas
- Bevelled edges for trip-proof access
- Material: polyurethane (PU)
- Material thickness: 13 mm
- Suitable for all floor types

Light grey

Item no. 42840 Total size (two parts): approx. 180 x 64 cm



With anti-slip surface.

POSITIONING MATERIAL



POSITIONING MATERIAL

Hand-sewn top quality produced by us.



MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

Colour chart, see cover



Full rolls

- Made from high-quality foam shaped sections with high density
- Shape-consistent, tested compression strength
- With zip for convenient removal of the covers

Half-rolls

- Made from high-quality foam shaped sections with high density
- Shape-consistent, tested compression strength
- With zip for convenient removal of the covers

50 x 20 cm
Item no. 48811

60 x 15 cm
Item no. 40197

TOP SELLER

60 x 15 x 7.5 cm
Item no. 40194

60 x 18 x 9 cm
Item no. 40195

TOP SELLER


50 x 15 cm
Item no. 48810

60 x 20 cm
Item no. 40198

60 x 25 x 12.5 cm
Item no. 40196

50 x 25 x 12.5 cm
Item no. 48817




 MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

 MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT



MADE
IN
GERMANY

MADE
IN
GERMANY



75 x 50 x 30 cm
With two leg recesses


Item no. 48826

TOP
SELLER

Colour chart, see cover


Lymphatic drainage wedges

- Made from high-quality foam shaped sections with high density
- Shape-consistent, tested compression strength
- With zip for convenient removal of the covers




60 x 50 x 20 cm
Without leg recess

Item no. 40877



75 x 50 x 30 cm
Without leg recess


Item no. 40878



75 x 22 x 30 cm
With one leg recess

Item no. 48825

TOP
SELLER



60 x 50 x 20 cm
With two leg recesses

Item no. 40872



75 x 55 x 20 cm
Without leg recess

Item no. 48821

Colour chart, see cover

Leg support cushion

- Made from high-quality foam shaped sections with high density
- Shape-consistent, tested compression strength
- With zip for convenient removal of the covers



SOCIAL RESPONSIBILITY

Collaboration with inclusive businesses to strengthen society.

For our bed production in Büren, we work with the Caritas Werkstätten St. Nikolaus. People with disabilities make metal and plastic production parts for us, prepare our practice stools for shipping, cut protective mats for transporting the therapy beds and lend a hand in many other areas.





MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

MADE IN GERMANY



Lordosis support cushions

- Made from high-quality foam shaped sections with high density
- Shape-consistent, tested compression strength



Massage cushions

- with soft foam filling

50 x 21.5 x 6.5 cm
Cervical spine
Item no. 40199

30 x 28 x 5 cm
Lumbar spine
Item no. 48861 **TOP SELLER**

ø 32 cm
Item no. 48832

40 x 40 cm
Item no. 48835



45 x 20 x 9 cm
Item no. 48805

40 x 30 x 9 cm
Item no. 48833



MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

MADE IN GERMANY



Crescent-shaped cushion

- Made from polystyrene core with foam cover
- High resistance to twisting



Positioning cube

- Made from polystyrene core with foam cover
- High resistance to twisting

80 x 40 x 40 cm
Asymmetrical
Item no. 48820



50 x 45 x 40 cm
Item no. 48801

50 x 25 x 20 cm
Item no. 48800



60 x 40 x 30 cm
Item no. 48802 **TOP SELLER**

60 x 50 x 40 cm
Item no. 48803

THE MULTITALENT

Ideal for patient positioning and fixation.

PINOFIT® Triblox Dark grey set of 3

- Set of three from shape-consistent, high-quality and easy to clean polyurethane
- For the easy and comfortable positioning of extremities, back and cervical spine
- Suitable for the fixation of joints during manual therapies
- Ideal for the positioning of feet, knees and arms

Item no. 43125
Each approx. 22 cm long, 15 cm wide and each approx. 7 to 9 cm tall



Arm positioning



Leg positioning



Foot positioning



Mobilisation wedges

- Made from polyurethane
- Stops slipping on the underlay
- For the treatment of extremities and spinal joints in the clinic and for patients' home exercises

Dark grey
Item no. 48830 Approx. 16 x 11 x 5 cm

Dark grey
Item no. 48831 Approx. 20 x 11 x 7 cm



Comfort full bolster extra soft with polyurethane cover

- Filled with extremely soft and adaptable polyester fibre fleece
- The positioning bolster perfectly fits body contours and allows highly comfortable treatment
- The cover made out of skin-friendly polyurethane is pleasant for the hands and can easily be wiped and disinfected
- Cover is washable at 60°C

Item no. 48834 Approx. 63 x 18 cm, cover white



PRONE POSITION CUSHION

For maximum lying comfort during treatment.



Anatomically moulded headrest

✓ MEDICAL PRODUCT CE

MADE OF ANTIMICROBIAL IMITATION LEATHER DISINFECTABLE WITH ALCOHOL-FREE DISINFECTANT

PINO prone position cushion

- Supports and facilitates therapy in prone position
- Nose and mouth are free; breathing is assisted by air slits at the front and on the sides
- The cushion can be adjusted to individual positioning requirements according to shoulder width and head size by varying the position of the bevelled upper body bolster
- Material from comfort foam with latex overlay and synthetic leather cover
- Easy to clean and simple to disinfect (e.g. with Septapin® Aquaquick p. 104)

All colours freely combinable

Lime
Item no. 40810 Approx. 46 x 35 x 10.5 cm, incl. cover

Dark grey
Item no. 40811 Approx. 46 x 35 x 10.5 cm, incl. cover



Anti-slip underside.


GEL CUSHION

Maximum comfort for your patients.



PINO gel pillow pearls, 3-chamber

- Particularly ergonomic thanks to its shape-consistent gel pearl filling
- Can be adjusted to any facial contour thanks to an innovative 3-chamber system, thus preventing uncomfortable pressure marks
- The separation of the cushion into one forehead and two cheek chambers helps prevent the shifting of the special gel filling
- Ideally combined with disposable nose slit tissues for hygienic comfort
- Easy to clean with soap and water or with the alcohol-free SEPTAPIN® AQUAQUICK rapid disinfection p. 104



New

Our highlights for your application

- ⊕ With 3-chamber system for even more comfort: 1. forehead chamber + 2. left cheek + 3. right cheek
- ⊕ Face lies even more evenly
- ⊕ Ideally combined with disposable nose slit tissues for hygienic comfort

Orange red

Item no. 40834 Approx. ø 30 cm

Lime azure

Item no. 40837 Approx. ø 30 cm

Water blue

Item no. 40838 Approx. ø 30 cm

Transparent

Item no. 40830 Approx. ø 30 cm

Saffron

Item no. 40839 Approx. ø 30 cm

Pink purple

Item no. 40841 Approx. ø 30 cm

PINOPRO gel cushions Pro

- Particularly ergonomic thanks to special gel filling
- Can be adjusted to any facial contour thanks to an innovative 2-chamber system, thus preventing uncomfortable pressure marks
- The separation of the cushions into separate forehead and cheek chambers helps prevent the shifting of the special gel filling
- Ideally combined with disposable nose slit tissues for hygienic comfort

Transparent

Item no. 40821 Approx. ø 30 cm



NOSE SLIT PILLOW

With gel inlay for optimal support.



With gel inlay

Nose slit pillow with gel inlay

- Comfortable face support
- Made from viscose elastic foam and gel inlay
- Cover from easy care, soft polyester
- Washable at 30°C

Grey

Item no. 40813 Approx. 30 x 26 x 5 cm

Lime

Item no. 40812 Approx. 30 x 26 x 5 cm



Spare cover for nose slit pillow with gel inlay

- Cover optionally made of easy-care, soft polyester or soft water-repellent polyurethane
- Washable at 30°C
- Concealed zip

Grey, polyester

Item no. 40814 Approx. 30 x 26 x 5 cm

Lime, polyester

Item no. 40816 Approx. 30 x 26 x 5 cm

Grey, polyurethane

Item no. 40818 Approx. 30 x 26 x 5 cm

Lime, polyurethane

Item no. 40819 Approx. 30 x 26 x 5 cm



Padding colours for therapy beds, stools & positioning material

All artificial leather covers are antimicrobial and feature a different surface structure.¹



Please choose a standard frame colour



Colour chart for terry fitted sheets¹



¹ Please note that it is often not possible to replicate a colour exactly in print. It is therefore possible that colours are not correctly displayed in our print media. Therefore all colours are supplied without liability, and deviations are possible.